



# Surgical Instruments Catalogue



HOMEDICA  
INTERNATIONAL



## **HOMEDICA INTERNATIONAL**

Homedica feels great honor introducing itself to be one of the leading manufacturing company of surgical & dental instruments. Professional always prefer Homedica instruments, which is our pride.

Our all instruments are manufactured by the high quality stainless steel made on our requirements, every instrument is hand crafted with the highest workmanship available in the industry. Through heat treatment, the instruments are hardened and tempered for proper functioning and long life. All instruments are polished electrolytically and manually for maximum resistance to corrosion and long life. Homedica pays special attention to the needs of doctors, with innovative designs to meet the needs of professionals. The wide assortment of standard patterns and exclusive Homedica designs gives you the most desirable range of surgical and dental instruments.

Over the several years, we have earned the reputation of for maintaining highest quality standards of quality and service. To manufacture instruments that meet with your changing demands is our achievement. Our highly skilled staff of professionals will continue to serve your needs and will try their outmost to uplift the standards which we have set in the past.

With your co-operation and guidance, Homedica will continue serving surgical and dental instruments in future.

**Management of Homedica international**

## BASIC METHODS FOR HANDLING AND CARE OF SURGICAL AND DENTAL INSTRUMENTS

Dear Customers:

We have received some questions concerning staining, pitting or rusting in stainless steel instruments. So, we took great honour and pleasure in guiding our respected customers to eliminate the problems about staining or rusting etc. The points are as follows:-

- 1) Rinse the instruments as soon as possible after use to remove blood and serum.
- 2) Dry the instruments after cleaning or before wrapping and storing, if not so, moisture left can cause rust and corrosion.
- 3) Do not leave instruments in cold sterilization solution for a long period, as the solution may lead the instrument to be corrosive.
- 4) Never use caustic cleaners, use those specified for cleaning surgical instruments.
- 5) Avoid cleaning instruments with steel wool; if want to clean with hand, use detergent soap and a soft brush and finally treat with distilled water.
- 6) You should check your autoclave; Iron, calcium, sodium or copper in water can cause rusting, staining or corrosion.
- 7) Leaving of instruments for a long time in the following solutions can also cause corrosion.
  - i) Phenol or Iodine.
  - ii) Aluminium, Barium or Calcium.
  - iii) Any Acid, Mercury or Potassium solution.
- 8) Try your utmost not to sterilize stainless steel instrument with chrome plated brass instruments, as any break in chrome surface may cause pitting or rusty surface in instruments.
- 9) Muslin wraps should be rinsed sufficiently by the laundry, so that all chlorine bleach may remove, which can cause corrosion or custing during autoclaving.

We hope that by acting on the above instructions, your problems relating to rusting, staining or pitting would abolish and the age of your instrument would increase.

**MANAGEMENT**

# ENT Surgery





HI-10-001  
4.5 mm Ø



HI-10-002  
5.5 mm Ø



HI-10-003  
6.5 mm Ø



HI-10-004  
7.5 mm Ø



HI-10-007  
2 mm Ø



HI-10-008  
3 mm Ø



HI-10-009  
4 mm Ø



HI-10-010  
5 mm Ø



**Toynbee Ear Speculum**  
3.6 cm 1 1/2"  
Set for adults

HI-10-005  
Set of 3  
5.5 mm  
6.5 mm  
7.5 mm

HI-10-006  
Set of 4  
4.5 mm  
5.5 mm  
6.5 mm  
7.5 mm



**Poltzer Ear Speculum**  
Set for adults

HI-10-011  
Set of 3  
2 mm  
3 mm  
4 mm

HI-10-012  
Set of 4  
2 mm  
3 mm  
4 mm  
5 mm



**Gruber Ear Speculum**  
Set for adults



HI-10-013  
a - 4.5 mm Ø



HI-10-014  
a - 5.5 mm Ø



HI-10-015  
a - 6.5 mm Ø



HI-10-016  
a - 7.5 mm Ø



**Gruber Ear Speculum**  
Set for adults



HI-10-017  
a - 4.5 mm Ø



HI-10-018  
a - 5.5 mm Ø



HI-10-019  
a - 6.5 mm Ø



HI-10-020  
a - 7.5 mm Ø



HI-10-021  
4.5 mm Ø



HI-10-022  
5.5 mm Ø



HI-10-023  
6.5 mm Ø



HI-10-024  
7.5 mm Ø



HI-10-025  
3 mm Ø



HI-10-026  
4.5 mm Ø



HI-10-027  
5.5 mm Ø



**Boucheron Ear Speculum**  
Set for adults



**Boucheron Ear Speculum**  
Set for children

## TUNING FORKS

Stimmgabeln | Diapasones  
Diapasons | Diapasoni



Hartmann  
Tuning Fork  
HI-10-028  
C 64



Hartmann  
Tuning Fork  
HI-10-029  
c 128



Hartmann  
Tuning Fork  
HI-10-030  
c1 256



Hartmann tuning  
fork with or without  
plastic base a1 440  
HI-10-031  
a1 440




Tuning Forks neurolog  
with hammer  
HI-10-032

## NOISE DRUM

Lärmtrommel | Aparato  
Assourdisseur | Assordatore



Barany Noise Drum  
HI-10-033



Lucae  
Tuning Fork  
HI-10-034  
c 128  
Variable up to h (ehn, b)



Rydel-Seiffer  
Tuning Fork  
HI-10-035  
C 64 / c 128  
With clamps C 64,  
Without clamps c 128



French  
Tuning Fork  
HI-10-036  
C2 16



French  
Tuning Fork  
HI-10-037  
C1 32



HI-10-038  
c1 128

HI-10-039  
c1 256

HI-10-040  
c2 512

HI-10-041  
c3 1024

HI-10-042  
c4 2048

HI-10-043  
c5 4096

Hartmann Tuning Fork



HI-10-044  
Hartmann tuning Forks Set case/5



HI-10-045  
c1 128

HI-10-046  
c1 256

HI-10-047  
c2 512

HI-10-048  
c3 1024

HI-10-049  
c4 2048



HI-10-050  
TUNING FORKS Set in case of 5

# EAR CURETTES



Ohrküretten

Curetas auriculares

Curettes auriculaires

Curette auricolari

	<b>HI-10-051</b> 16.5 cm 6 1/2" Blunt	<b>HI-10-052</b> 16.5 cm 6 1/2" Blunt	<b>HI-10-053</b> 16.5 cm 6 1/2" Blunt	<b>HI-10-054</b> 16.5 cm 6 1/2" Blunt	<b>HI-10-055</b> 16.5 cm 6 1/2" Blunt	<b>HI-10-056</b> 16.5 cm 6 1/2" Blunt	<b>HI-10-057</b> 16.5 cm 6 1/2" Blunt	<b>HI-10-058</b> 16.5 cm 6 1/2" Blunt
	<b>HI-10-059</b> 16.5 cm 6 1/2" Sharp	<b>HI-10-060</b> 16.5 cm 6 1/2" Sharp	<b>HI-10-061</b> 16.5 cm 6 1/2" Sharp	<b>HI-10-062</b> 16.5 cm 6 1/2" Sharp	<b>HI-10-063</b> 16.5 cm 6 1/2" Sharp	<b>HI-10-064</b> 16.5 cm 6 1/2" Sharp	<b>HI-10-065</b> 16.5 cm 6 1/2" Sharp	<b>HI-10-066</b> 16.5 cm 6 1/2" Sharp
<b>HI-10-067</b> 16.5 cm 6 1/2" Curved, Blunt	<b>HI-10-068</b> 16.5 cm 6 1/2" Curved, Blunt	<b>HI-10-069</b> 16.5 cm 6 1/2" Curved, Blunt	<b>HI-10-070</b> 16.5 cm 6 1/2" Curved, Blunt	<b>HI-10-071</b> 16.5 cm 6 1/2" Curved, Blunt	<b>HI-10-072</b> 16.5 cm 6 1/2" Curved, Blunt	<b>HI-10-073</b> 16.5 cm 6 1/2" Curved, Blunt	<b>HI-10-074</b> 16.5 cm 6 1/2" Curved, Blunt	
<b>HI-10-075</b> 16.5 cm 6 1/2" Curved, Sharp	<b>HI-10-076</b> 16.5 cm 6 1/2" Curved, Sharp	<b>HI-10-077</b> 16.5 cm 6 1/2" Curved, Sharp	<b>HI-10-078</b> 16.5 cm 6 1/2" Curved, Sharp	<b>HI-10-079</b> 16.5 cm 6 1/2" Curved, Sharp	<b>HI-10-080</b> 16.5 cm 6 1/2" Curved, Sharp	<b>HI-10-081</b> 16.5 cm 6 1/2" Curved, Sharp	<b>HI-10-082</b> 16.5 cm 6 1/2" Curved, Sharp	

Buck, Ear  
Curette



Pictures show the outer diameter





**Billeau  
Ear Loop**



**HI-10-083**  
16.5 cm 6 1/2"



**HI-10-084**  
16.5 cm 6 1/2"



**HI-10-085**  
16.5 cm 6 1/2"



**Lucae  
Hooklet**



**HI-10-086**  
14 cm 5 1/2"



**HI-10-087**  
14 cm 5 1/2"



**HI-10-088**  
14 cm 5 1/2"



**Wagener  
Hooklet Buttoned**



**HI-10-089**  
14 cm 5 1/2"



**HI-10-090**  
14 cm 5 1/2"



**HI-10-091**  
14 cm 5 1/2"



**HI-10-092**  
14 cm 5 1/2"



**HI-10-093**  
14 cm 5 1/2"



**HI-10-094**  
14 cm 5 1/2"



**HI-10-095**  
14 cm 5 1/2"



**Farrell Applicator**  
**HI-10-096**  
16.5 cm 6 1/2"  
cross-serrated tip

# EAR POLYPUS SNARES



*Ohrpolypenschnürer*  
*Serre noeuds auriculaires*

*Polipótomos auriculares*  
*Serranodi auricolari*



**Nasal Polypus Snares**  
**HI-10-097**  
0.3 mm Ø X 100 mm  
In bundles of 12 pieces



**Krause Nasal Polypous Snares**  
**HI-10-098**  
25.4 cm 10"



**Kraus -Vosse Nasal Polypous Snares**  
**HI-10-099**  
25.4 cm 10"  
With 2 inserts



**Kraus -Wilde Nasal Polypous Snares**  
**HI-10-100**  
24 cm 9 1/2"



**HI-10-101**  
50CC Complete  
**HI-10-102**  
75CC Complete

**HI-10-103**  
100CC Complete  
**HI-10-104**  
150CC Complete  
**Reiner Ear syringe**  
With Tip & Shield



**Reiner Ear Syringe Tip Only**  
**HI-10-105**  
Tip Only



**Reiner Ear Syringe Tip Only**  
**HI-10-106**  
Tip Only



**Reiner Ear Syringe Shield**  
**HI-10-107**



**Littauer Ear  
Polypus Forceps**  
HI-10-108  
12 cm 4 3/4"



**Hartmann Ear  
Polypus Forceps**  
HI-10-109  
12 cm 4 3/4"  
Delicate jaw



**Hartmann Ear  
Polypus Forceps**  
HI-10-110  
14 cm 5 1/2"  
Normal jaw



**Quire Ear  
Polypus Forceps**  
HI-10-111  
10.5 cm 4 1/4"  
Foreign body lever



**Buck Ear  
Polypus Forceps**  
HI-10-112  
11.5 cm 4 1/2"  
Foreign body forceps

# EAR FORCEPS



Ohrzangchen | Pinzas auriculares  
Pincers auriculaires | Pinze auricolari



**Hartman  
Ear Forceps**  
HI-10-113



**Hartmann  
Ear Forceps**  
HI-10-114  
X = 8.5 cm, 3 1/4"  
HI-10-115  
X = 14 cm, 5 1/2"



**Hartmann  
Ear Forceps**  
HI-10-116  
X = 8.5 cm, 3 1/4"  
HI-10-117  
X = 14 cm, 5 1/2"  
HI-10-118  
X = 16 cm, 6 1/4"



**Hartmann  
Ear Forceps**  
HI-10-119 2 mm Ø  
X = 8.5 cm, 3 1/4"  
HI-10-120 2.5 mm Ø  
X = 8.5 cm, 3 1/4"  
HI-10-121 3 mm Ø  
X = 8.5 cm, 3 1/4"



**Hartmann  
Ear Forceps**  
HI-10-122



**Hartmann-Wullstein  
Ear Forceps**  
HI-10-123  
X = 8.5 cm, 3 1/4"



**Hartmann-Wullstein  
Ear Forceps**  
HI-10-124  
1 mm Ø  
X = 8.5 cm, 3 1/4"



**Hartmann-Weingartner  
Ear Forceps**  
HI-10-125  
X = 8.5 cm, 3 1/4"



**Duckbill  
Ear Forceps**  
HI-10-126



**Weingartner  
Ear Forceps**  
HI-10-127



**Strumpel  
Ear Forceps**  
HI-10-128



**Hartmann-hoffmann  
Ear Forceps**  
HI-10-129



**Hartmann-Noyes  
Ear Forceps**  
HI-10-130  
X = 9 cm, 3 1/2"



**Struempel  
Ear Forceps**  
HI-10-131 0°  
X = 8.5 cm, 3 1/4"



**Struempel  
Ear Forceps**  
HI-10-132 45°  
X = 8.5 cm, 3 1/4"




**Struempel  
Ear Forceps**  
HI-10-133 90°  
X = 8.5 cm, 3 1/4"



**MICRO CUP FORCEPS**  
HI-10-134  
X = 75 mm 3"  
a = 4.0 mm  
b = 1.0 mm  
c = 0.9 mm  
Straight, Ovel

**MICRO CUP FORCEPS**  
HI-10-135  
X = 75 mm 3"  
a = 4.0 mm  
b = 1.0 mm  
c = 0.9 mm  
Curved upwards, Ovel

**MICRO CUP FORCEPS**  
HI-10-136  
X = 75 mm 3"  
a = 4.0 mm  
b = 1.0 mm  
c = 0.9 mm  
Curved to left, Ovel



**MICRO CUP FORCEPS**

**MICRO CUP FORCEPS**  
HI-10-137  
X = 75 mm 3"  
a = 4.0 mm  
b = 1.0 mm  
c = 0.9 mm  
Curved to right, Ovel


**MICRO CUP FORCEPS**  
HI-10-138  
X = 75 mm 3"  
a = 4.0 mm  
b = 1.0 mm  
c = 0.9 mm  
Curved downwards, Ovel



**Citelli Bone  
 Ear Rongeurs**  
 HI-10-139  
 17 cm 6 3/4"  
 X = 60 mm



**Kerrison Bone  
 Ear Rongeurs**  
 HI-10-140  
 20 cm 8"  
 X = 80 mm



**Beyer Bone  
 Ear Rongeurs**  
 HI-10-141  
 14 cm 5 1/2"  
 1 mm



**Beyer Bone  
 Ear Rongeurs**  
 HI-10-142  
 19 cm 7 1/2"  
 2 mm

# NASAL SPECULA



*Nasenspekula* | *Especulos nasales*  
*Spéculums nasaux* | *Specoli nasali*



**Thudichum  
(Goldsmith)  
Nasal Speculum**



**HI-10-143**  
6.5 cm 2 1/2"  
21 mm X 9.5 mm




**HI-10-144**  
6.5 cm 2 1/2"  
21 mm X 10.5 mm



**HI-10-145**  
6.5 cm 2 1/2"  
21 mm X 11.5 mm




**Hartmann  
Nasal Speculum**



**HI-10-146**  
16 cm 6 1/4"



**HI-10-147**  
16 cm 6 1/4"



**HI-10-148**  
16 cm 6 1/4"



**Hartmann-Halle  
Nasal Speculum**



**HI-10-149**  
15 cm 6"



**HI-10-150**  
15 cm 6"



**HI-10-151**  
15 cm 6"

# NASAL SPECULA

Nasenspekula

Especulos nasales

Spéculums nasaux

Specoli nasali



**HI-10-152**  
13.5 cm 5 1/4"  
For babies

**Tieck-Halle**  
Nasal Speculum




**HI-10-153**  
14 cm 5 1/2"


**Schwarz Nasal Speculum**



**Vienna Speculum**  
**HI-10-154**  
14.6cm 5 3/4  
medium standard patter




**HI-10-155**  
13 cm 5"  
X = 35 mm




**HI-10-156**  
13 cm 5"  
X = 50 mm



**Killian Nasal Speculum**



**HI-10-157**  
13 cm 5"  
X = 75 mm



**HI-10-158**  
13 cm 5"  
X = 90 mm

# NASAL SPECULA



*Nasenspekula* | *Especulos nasales*  
*Spéculums nasaux* | *Specoli nasali*



**Killian Septum Speculum**  
**HI-10-159**  
 15.2 cm 6"  
 2.51 cm blades  
 withset screw




**Killian Septum Speculum**  
**HI-10-160**  
 15.2 cm 6"  
 8.9cm 3 1/2 blades  
 set screw




**Killian Septum Speculum**  
**HI-10-161**  
 Fiber Optic 5 1/2" 14.0mm  
 length 2.51mm deep




**Beckmann**  
**Nasal Speculum**



**HI-10-162**  
 15 cm 6"




**HI-10-163**  
 15 cm 6"




**HI-10-164**  
 15 cm 6"




**Mod. Wien**  
**Nasal Speculum**



**HI-10-165**  
 14.5 cm 5 3/4"



**HI-10-166**  
 14.5 cm 5 3/4"



**HI-10-167**  
 14.5 cm 5 3/4"



## NASAL SPECULA

Nasenspekula


Especulos nasales

Spéculums nasaux


Specoli nasali




Cottle Nasal Speculum




**HI-10-168**  
15 cm 6"  
X = 35 mm



**HI-10-169**  
15 cm 6"  
X = 50 mm




**HI-10-170**  
15 cm 6"  
X = 75 mm




**HI-10-171**  
15 cm 6"  
X = 90 mm



Cottle Nasal Speculum



**HI-10-172**  
14 cm 5 1/2"  
X = 45 mm



**HI-10-173**  
14 cm 5 1/2"  
X = 75 mm

## NASAL POLYPUS AND SEPTUM FORCEPS



Nasenpolypenzangen-Septumzangen  
Pincas à polypes et la cloison nasale

Pinzas para pólipos y el tabique nasal  
Pinze per polipi e setto nasale



**Westmacott Nasal Polypus and Septum Forceps**  
HI-10-174 20 cm 8"



**Lubet-Barbon Nasal Polypus and Septum Forceps**  
HI-10-175 20 cm 8"



**Hartmann Nasal Polypus and Septum Forceps**  
HI-10-176 16 cm 6 1/4"  
HI-10-177 20 cm 8"



**Tilley Nasal Polypus and Septum Forceps**  
HI-10-178 14 cm 5 1/2"



**Heymann-Knight Nasal Polypus and Septum Forceps**  
HI-10-179 18 cm 7"



**Craig Nasal Polypus and Septum Forceps**  
HI-10-180 18 cm 7"



**Knight Nasal Polypus and Septum Forceps**  
HI-10-181 18 cm 7"



**Littauer Nasal Polypus and Septum Forceps**  
HI-10-182 18 cm 7"

Nasenpolypenzangen, Nasentamponzange


Pinzas para pólipos nasales, Pinza nasale para tamponiamento

Pinces à polypes nasales, Pince à pansement pour le nez

Pinze per polipi nasali, Pinza per tamponamento nasale



**Heymann Nasal Dressing Forceps**



**Bruenings Nasal Cutting Forceps**




**Takahashi Nasal Cutting Forceps**



**Noyes Nasal Dressing Forceps**

**HI-10-192**  
 X = 11 cm 3 1/2"  
 1 X 2



**Luc Nasal Cutting Forceps**



**Weil Nasal Cutting Forceps**

**Watson-Williams Nasal Cutting Forceps**

**HI-10-199**  
 X = 12 cm 4 3/4"

# TONSIL KNIVES



Tonsillennmesser

Bisturis para amigdales

Couteaux à tonsils

Bisturi per tonsille



**Freer  
Tonsil Knife  
HI-10-200**  
21 cm 8 1/4"



**Freer  
Tonsil Knife  
HI-10-201**  
21 cm 8 1/4"



**Tonsil Knife  
HI-10-202**  
20 cm 8"  
Sharp



**Tonsil Knife  
HI-10-203**  
20 cm 8"  
Sickle knife



**Bruenings  
Tonsil Knife  
HI-10-204**  
23 cm 9"



**Abraham  
Tonsil Knife  
HI-10-205**  
22 cm 8 3/4"



**Fischer Tonsil  
Knife and Dissector  
HI-10-206**  
21.6 cm 8 1/2"  
semi-sharp tip, with teeth



**EVE Tonsil Snare**  
27.9 cm 11"  
**HI-10-207** without ratchet  
**HI-10-208** with ratchet

Septumelevatorien - Raspatorien

Elevadores - legras para el tabique nasal

Elevateurs - rugines pour la cloison nasale

Elevatori - periostotomi per setto



**Lempert Septum Elevators Raspatory HI-10-209**  
20.5 cm 8"  
X = 6 mm



**Lempert Septum Elevators Raspatory HI-10-210**  
18 cm 7"  
X = 6 mm



**Cottle Septum Elevators Raspatory HI-10-211**  
19 cm 7 1/2"  
X = 8 mm  
Sharp



**Joseph Septum Elevators Raspatory HI-10-212**  
16 cm 6 1/4"  
X = 5 mm  
Periosteal elevators



**Joseph-Masing Septum Elevators Raspatory HI-10-213**  
21.5 cm 8 1/2"  
Sharp



**Pennington Septum Elevators Raspatory HI-10-214**  
21.5 cm 8 1/2"  
Blunt



**Freer Septum Elevators Raspatory HI-10-215**  
19 cm 7 1/2"  
Blunt



**Freer Septum Elevators Raspatory HI-10-216**  
18 cm 7"  
Sharp/Blunt



**Freer Septum Elevators Raspatory HI-10-217**  
20.5 cm 8"  
Sharp/Blunt



**Cottle Septum Elevators Raspatory HI-10-218**  
23.5 cm 9 1/4"  
X = 5 mm  
Sharp

# PERIOSTEAL ELEVATORS, SEPTUM ELEVATORS



Periostelevatorien, Septumelevatorien

Legras para periostio, Elevadores para el tabique nasal

Elevateurs pour la périoste, Elevateurs pour la cloison nasale

Elevatori per il periostio, Elevatori per setto



**Masing Periosteal Septum Elevators**  
HI-10-219  
20 cm 8"




**Ballenger-Hajek Periosteal Septum Elevators**  
HI-10-220  
19 cm 7 1/4"  
HI-10-221  
22 cm 8 3/4"



**Cottle Periosteal Septum Elevators**  
HI-10-222  
22 cm 8 3/4"



**Cottle Periosteal Septum Elevators**  
HI-10-223  
21 cm 8 1/4"



**Watson-Cheyne**  
HI-10-224  
13 cm 5"  
Raspatory



**Boies**  
HI-10-225  
19 cm 7 1/2"  
Nasal fracture elevator,  
Blunt



**Goldman**  
HI-10-226  
17 cm 6 3/4"  
Septum displacer

## NASAL SAWS

*Nasensägen*

*Sirerras nasales*

*Scies nasales*

*Seghe nasali*



**Joseph Nasal Saw**  
HI-10-227  
19 cm 7 1/2"



**Joseph Nasal Saw**  
HI-10-228  
19 cm 7 1/2"



**Joseph Nasal Saw**  
HI-10-229  
19 cm 7 1/2"



**Maltz Rasp**  
HI-10-230  
21 cm 8-1/4"  
improved double end  
pattern, forward and  
backward cutting



**Joseph Rasp**  
HI-10-231  
15.9 cm 6-1/4"  
fine cross serrations



**Joseph Nasal Saw**  
HI-10-232  
17 cm 6 3/4"

**Joseph Nasal Saw**  
HI-10-233  
17 cm 6 3/4"

# RASPS WITH TUNGSTEN CARBIDE



Raspeln mit Hartmetall

limas con carburo de tungsteno

Rapes avec inserts de carbure de tungstene

raspe con tungsteno




**Fomon  
Rasp with  
Tungsten  
Carbide**



**HI-10-234**  
20 cm 8"  
Teeth 12




**HI-10-235**  
20 cm 8"  
Teeth 16




**HI-10-236**  
20 cm 8"  
Teeth 24



**HI-10-237**  
20 cm 8"  
Teeth 36



**HI-10-238**  
20 cm 8"  
Teeth 52



**Fomon Rasp  
with Tungsten  
Carbide TC**



**HI-10-239**  
20 cm 8"  
Teeth 14



**HI-10-240**  
20 cm 8"  
Teeth 20



**HI-10-241**  
20 cm 8"  
Teeth 30



**HI-10-242**  
20 cm 8"  
Teeth 44




**HI-10-243**  
20 cm 8"  
Teeth 62




# NASAL SAWS

Nasensägen | Sirerras nasales  
Scies nasales | Seghe nasali



**Cottle Nasal Saw**  
HI-10-244  
21 cm 8 1/4"

**Maltz Nasal Saw**



**HI-10-245**  
17 cm 6 3/4"  
34 mm x 7 mm  
Pushing cut



**HI-10-246**  
17 cm 6 3/4"  
34 mm x 7 mm  
Drawing cut



**Aufricht Nasal Saw**  
Glabella Rasp,  
backward cutting



**HI-10-247**  
21 cm 8 1/4"  
Pushing cut




**HI-10-248**  
21 cm 8 1/4"  
Drawing cut



**HI-10-249**  
17.5 cm 7"  
40 mm X 9 mm



**Wiener Nasal Saw**



**HI-10-250**  
18 cm 7"  
3.5 mm Ø  
Antrum rasps



**HI-10-251**  
18 cm 7"  
3.5 mm Ø  
Antrum rasps



**Gallaher Nasal Saw**



**HI-10-252**  
16 cm 6 1/4"



**HI-10-253**  
16 cm 6 1/4"



**HI-10-254**  
16 cm 6 1/4"



**Fomon Nasal Saw**  
HI-10-255  
21 cm 8 1/4"



# SEPTUM INSTRUMENTS



Septuminstrumente  
instruments pour la cloison nasal

Instrumentos para el tabique nasal  
Strumenti per il setto nasale



**Bruenings-Brunk**  
**HI-10-256**  
19 cm 7 1/2"  
Antrum Retractor  
Serrated



**Ritter-Halle**  
**HI-10-257** 16 cm 6 1/4"  
2.5 mm Ø  
**HI-10-258** 16 cm 6 1/4"  
3.0 mm Ø  
**HI-10-259** 16 cm 6 1/4"  
4.0 mm Ø  
Fontal sinus bougies



**Cottle**  
**HI-10-260**  
1.2 cm 4 3/4"  
Handle



**Cottle**  
**HI-10-261**  
7 mm  
Guide support  
for osteotomes



**Cottle**  
**HI-10-262**  
65 mm X 30 mm X 30 mm  
Cartilage crusher



**Cottle-Walsham**  
**Septum Forceps**




**HI-10-263**  
23 cm 9"  
Straight



**HI-10-264**  
23 cm 9"  
Left



**HI-10-265**  
23 cm 9"  
Right



**Ash Septum**  
**Forceps**  
**HI-10-266**  
23 cm 9"

## SEPTUM FORCEPS

Septumzangen

Pinces pour la cloison nasale

Pinzas para el tabique nasal


Pinze per il setto nasale



**Middleton-Jansen  
Septum Forceps**  
HI-10-267  
19 cm 7 1/2"



**Middleton-Jansen  
Septum Forceps**  
HI-10-268  
19 cm 7 1/2"



**Caplan Septum  
Forceps**  
HI-10-269  
20 cm 8"  
Toothed



**Rubin Septum  
Forceps**  
HI-10-270  
21 cm 8 1/4"  
Protection cap

*adenoid currettes*

*Cuchillos para vegetaciones adenoides*

*Couteaux a vegetations adenoides*

*Adenotomi*



**HI-10-271**  
21 cm 8 1/4"



**HI-10-272**  
21 cm 8 1/4"



**HI-10-273**  
21 cm 8 1/4"



**HI-10-274**  
21 cm 8 1/4"



**HI-10-275**  
22 cm 8 3/4"



**HI-10-276**  
22 cm 8 3/4"



**HI-10-277**  
22 cm 8 3/4"



**HI-10-278**  
22 cm 8 3/4"

**St. Clair-Thomson  
Adenoid Curette**

**Barnhill  
Adenoid Curette**



**HI-10-279**  
22 cm 8 3/4"  
Straight



**HI-10-280**  
22 cm 8 3/4"  
Straight



**HI-10-281**  
22 cm 8 3/4"  
Straight



**HI-10-282**  
22 cm 8 3/4"  
Straight



**HI-10-283**  
22 cm 8 3/4"  
Straight



**HI-10-284**  
22 cm 8 3/4"  
Straight



**HI-10-285**  
22 cm 8 3/4"  
Curved



**HI-10-286**  
22 cm 8 3/4"  
Curved



**HI-10-287**  
22 cm 8 3/4"  
Curved



**HI-10-288**  
22 cm 8 3/4"  
Curved



**HI-10-289**  
22 cm 8 3/4"  
Curved




**HI-10-290**  
22 cm 8 3/4"  
Curved

**Beckmann  
Adenoid Curette**


**Beckman  
Adenoid Curette**




La Force




La Force  
adenotome  
unperforated  
Fig 1  
**HI-10-291**  
21.6 cm 8 1/2"  
25 x 16.0mm  
Infant




La Force  
adenotome  
perforated  
Fig 1  
**HI-10-292**  
21.6 cm 8 1/2"  
25 x 16.0mm  
Infant




La Force  
adenotome  
unperforated  
Fig 2  
**HI-10-293**  
24.1 cm 9 1/2"  
26 x 19.3mm  
Medium



La Force  
adenotome  
perforated  
Fig 2  
**HI-10-294**  
24.1 cm 9 1/2"  
26 x 19.3mm  
Medium



La Force  
adenotome  
unperforated  
Fig 3  
**HI-10-295**  
25 cm 10"  
30 x 22.2mm  
Large



La Force  
adenotome  
perforated  
Fig 3  
**HI-10-296**  
25 cm 10"  
30 x 22.2mm  
Large



**HI-10-297**  
20 cm 8"  
X = 3 mm  
Straight  
**HI-10-298**  
20 cm 8"  
Straight  
X = 4 mm  
**HI-10-299**  
20 cm 8"  
Straight  
X = 5 mm

Ballenger Knife



**HI-10-300**  
20 cm 8"  
X = 3 mm  
Bayonet  
**HI-10-301**  
20 cm 8"  
Bayonet  
X = 4 mm  
**HI-10-302**  
20 cm 8"  
Bayonet  
X = 5 mm

Ballenger Knife



Wieder Tongue Depressor  
**HI-10-303** large blade 36mm  
wide serrated

Wieder Tongue Depressor  
**HI-10-304** 43.5mm width blade

# MOUTH GAGS

Mundsperrer | Abrebocas  
Ouvre-bouches | Apribocca



**Whitehead  
Mouth Gag**  
HI-10-305 11 cm 4 1/4"  
HI-10-306 13 cm 5"  
HI-10-307 15 cm 6"



**Jennings  
Mouth Gag**  
HI-10-308 9 cm 3 1/2"  
HI-10-309 11 cm 4 1/4"  
HI-10-310 13 cm 5"  
HI-10-311 15 cm 6"



**Heister  
Mouth Gag**  
HI-10-312  
13 cm 5"



**Mauder  
Mouth Gag**  
HI-10-313  
10 cm 4"  
4 mm Ø  
Wood



**Molt Mouth Gag**  
HI-10-314  
11 cm 4 1/4"  
For children



**Molt Mouth Gag**  
HI-10-315  
14 cm 5 1/2"  
For adults



**Denhart  
Mouth Gag**  
HI-10-316  
12 cm 4 3/4"



**Fergusson-Ackland  
Mouth Gag**  
HI-10-317  
14.5 cm 5 3/4"



**Roser-Konig  
Mouth Gag**  
HI-10-318  
16 cm 6 1/4"  
HI-10-319  
19 cm 7 1/2"



**Jackson**  
**HI-10-320**  
22 cm 8 3/4"  
Laryngeal forceps



**Marschik**  
**HI-10-321**  
22 cm 8 3/4"  
Tonsil seizing forceps




**Blohmke**  
**Tonsil**  
**Forceps**  
**HI-10-322**  
20 cm 8"



**Colver Tonsil**  
**Forceps**  
**HI-10-323**  
19 cm 7 1/2"  
Straight  
**HI-10-324**  
19 cm 7 1/2"  
Curved



**White Tonsil**  
**Seizing Forceps**  
**HI-10-325**  
19.1 cm 7 1/2"  
**HI-10-326**  
23 cm 9"



**Museux**  
**Tonsil Holding**  
**Forceps**  
**HI-10-327** 20 cm 8"

## TONSIL INSTRUMENTS



*Tonsilleninstrumente*

*Instrumentos para amígdalas*

*Instruments pour amygdales*

*Strumenti per tonsille*



**Wilde Nasal  
Dressing Forceps**  
HI-10-328  
16.5cm 6 1/2"  
angular



**Troeltsch  
Dressing Forceps**  
HI-10-329  
11cm, 4 1/2"



**Troeltsch  
Dressing Forceps**  
HI-10-330  
13 cm, 5"




**Lucae Dressing  
Forceps**  
HI-10-331  
14 cm, 5 1/2"



**Jansen Dressing  
Forceps**  
HI-10-332  
16 cm, 6 1/4"





**Tyding Tonsil Snare**  
**HI-10-333**  
 With one straight tip



**Bruenings**  
**HI-10-334**  
 10 m 33 ft  
 0.3 mm Ø  
**HI-10-335**  
 10 m 33 ft  
 0.4 mm Ø



**Bruenings**  
**HI-10-336**  
 40 mm  
 Package of 100 pieces  
**HI-10-337**  
 50 mm  
 Package of 100 pieces



**HI-10-338**  
 28 cm 11"

**HI-10-339**  
 28 cm 11"

**Eves Tonsil Snare**



**HI-10-340**  
 28 cm 11"

**Eves Tonsil Snare**



**HI-10-341**  
 28 cm 11"

**Bruenings Tonsil Snare**

# LARYNGEAL MIRRORS



*Kehlkopfspiegel* | *Especjos laringeos*  
*Miroirs laryngiens* | *Specchi laringei*



- HI-10-342 8 mm Ø
- HI-10-343 10 mm Ø
- HI-10-344 12 mm Ø
- HI-10-345 14 mm Ø
- HI-10-346 16 mm Ø
- HI-10-347 18 mm Ø
- HI-10-348 20 mm Ø
- HI-10-349 22 mm Ø
- HI-10-350 24 mm Ø
- HI-10-351 26 mm Ø
- HI-10-352 28 mm Ø
- HI-10-353 30 mm Ø

Laryngeal mirrors



Laryngeal Mirrors  
 Handle  
 HI-10-354  
 10 cm 4"



- HI-10-355 8 mm Ø
- HI-10-356 10 mm Ø
- HI-10-357 12 mm Ø
- HI-10-358 14 mm Ø
- HI-10-359 16 mm Ø
- HI-10-360 18 mm Ø
- HI-10-361 20 mm Ø
- HI-10-362 22 mm Ø
- HI-10-363 24 mm Ø
- HI-10-364 26 mm Ø
- HI-10-365 28 mm Ø
- HI-10-366 30 mm Ø

Laryngeal Mirrors  
 with Handle



- HI-10-367 10 mm Ø
- HI-10-368 12 mm Ø
- HI-10-369 14 mm Ø
- HI-10-370 16 mm Ø
- HI-10-371 18 mm Ø

Michel Moveable  
 Mirror

Trachealdilatatoren, Trachealkanülen

Dilatateurs pour la tracheotomie, Canules pour la tracheotomie

Dilatadores para la traqueotomia, Cabykas oara ka traqyea

Dilatatori tracheali, Cabbyke tracheali



**Laborde Trachea Dilator 3 Prong**  
HI-10-372  
14 cm 5 1/2"



**Trousseau Trachea Dilator 2 Prong**  
HI-10-373  
14 cm 5 1/2"



**Uckermann Trachea Tubes**  
HI-10-374  
13 cm 5"  
10 mm Ø



**Luer-Standard Trachea Tubes**  
German silver

- HI-10-375 4 mm Ø
- HI-10-376 5 mm Ø
- HI-10-377 6 mm Ø
- HI-10-378 7 mm Ø
- HI-10-379 8 mm Ø
- HI-10-380 9 mm Ø
- HI-10-381 10 mm Ø
- HI-10-382 11 mm Ø
- HI-10-383 12 mm Ø
- HI-10-384 13 mm Ø



**Chevalier-Jackson Trachea Tubes**  
German silver

- HI-10-385 4 mm Ø
- HI-10-386 5 mm Ø
- HI-10-387 6 mm Ø
- HI-10-388 7 mm Ø
- HI-10-389 8 mm Ø
- HI-10-390 9 mm Ø
- HI-10-391 10 mm Ø
- HI-10-392 11 mm Ø
- HI-10-393 12 mm Ø
- HI-10-394 13 mm Ø

# LARYNGEAL POLYPUD FORCEPS



*Kehlkopfpolypenzangen*  
*Pincas a polypes laryngiens*

*Pinzas para polipos laringeos*  
*Pinze per polipi laringei*

**Jako-Cherry Micro  
Laryngoscope**  
**HI-10-395**  
Type II, Newborn  
Small, Overall length 114mm  
min. I.D. 7.8mm



**Carlens Mediastinoscope, Short**  
**HI-10-396**  
Overall length 116mm, min. I.D. 10.8mm,  
nominal proximal dims 16.0mm, max O.D. 18.4mm



**Carlens Mediastinoscope, Long**  
**HI-10-397**  
Overall length 162mm, min. I.D. 11.0mm,  
nominal proximal dims 18.0mm, max O.D. 22.3mm



**Carlens Mediastinoscope,  
Extra-long**  
**HI-10-398**  
Overall length 191mm, min. I.D. 11.9mm,  
nominal proximal dims 19.0mm, max O.D. 22.3mm



**Tucker Mid-Lighted Side Slide  
Laryngoscop**  
**HI-10-399**  
Overall length 91mm, min. I.D. 6.9mm




**Tucker Mid-Lighted Side Slide  
Laryngoscope, Adolescent**  
**HI-10-400**  
Overall length 106mm, min. I.D. 9.0mm

**Tucker Mid-Lighted Side Slide  
Laryngoscope**  
**HI-10-401 Adult**  
Overall length 121mm, min. I.D. 10.7mm



**Tucker Mid-Lighted Side Slide  
Laryngoscope**  
**HI-10-402 Adult Long**  
Overall length 160mm, min. I.D. 11.8mm

(Note) Required size available on Demand



**Laryngoscope Holder HI-10-403**  
Chest Support complete, for adults and children

**Laryngoscope Holder HI-10-404**  
Chest Support complete, for adults and children

**Holinger Bronchoscope HI-10-405 3X20**  
Infant, Scope size 3 x 20, overall length 203mm, min. I.D. 2.9mm

**Holinger Bronchoscope HI-10-406 3.5X20**  
Infant, Scope size 3.5 x 20, overall length 203mm, min. I.D. 3.5mm

**Holinger Bronchoscope HI-10-407 9X40**  
Scope size 9 x 40, overall length 400mm, min. I.D. 8.5mm

**Window Plug For Bronchoscopes HI-10-408**  
1 1/8" 2.5cm

# LARYNGEAL POLYPUD FORCEPS



*Kehlkopfpolypenzangen*  
*Pinces a polypes laryngiens*

*Pinzas para polipos laringeos*  
*Pinze per polipi laringei*

**Jesburg Esophagoscope**  
**HI-10-409** 6 x 8 x32, Overall length 317mm, min. I.D. 5.1mm


**Jesburg Esophagoscope**  
**HI-10-410** 6.5 x 10 x 35, Overall length 349mm, min. I.D. 5.5mm

**Jesburg Esophagoscope**  
**HI-10-411** 7 x 10 x 42, Overall length 421mm, min. I.D. 5.8mm

**Suction Tube**  
**HI-10-412**  
with fix head

**Fenestrated Cupped Grasping Forceps**  
**HI-10-413**

**Frazier Suction Tube**  
**HI-10-414**  
11 French tip, Working length 19cm 7 1/2"



**HOMEDICA**  
INTERNATIONAL