


Surgical Instruments Catalogue



HOMEDICA
INTERNATIONAL



HOMEDICA INTERNATIONAL

Homedica feels great honor introducing itself to be one of the leading manufacturing company of surgical & dental instruments. Professional always prefer Homedica instruments, which is our pride.

Our all instruments are manufactured by the high quality stainless steel made on our requirements, every instrument is hand crafted with the highest workmanship available in the industry. Through heat treatment, the instruments are hardened and tempered for proper functioning and long life. All instruments are polished electrolytically and manually for maximum resistance to corrosion and long life. Homedica pays special attention to the needs of doctors, with innovative designs to meet the needs of professionals. The wide assortment of standard patterns and exclusive Homedica designs gives you the most desirable range of surgical and dental instruments.

Over the several years, we have earned the reputation of for maintaining highest quality standards of quality and service. To manufacture instruments that meet with your changing demands is our achievement. Our highly skilled staff of professionals will continue to serve your needs and will try their outmost to uplift the standards which we have set in the past.

With your co-operation and guidance, Homedica will continue serving surgical and dental instruments in future.

Management of Homedica international

BASIC METHODS FOR HANDLING AND CARE OF SURGICAL AND DENTAL INSTRUMENTS

Dear Customers:

We have received some questions concerning staining, pitting or rusting in stainless steel instruments. So, we took great honour and pleasure in guiding our respected customers to eliminate the problems about staining or rusting etc. The points are as follows:-

- 1) Rinse the instruments as soon as possible after use to remove blood and serum.
- 2) Dry the instruments after cleaning or before wrapping and storing, if not so, moisture left can cause rust and corrosion.
- 3) Do not leave instruments in cold sterilization solution for a long period, as the solution may lead the instrument to be corrosive.
- 4) Never use caustic cleaners, use those specified for cleaning surgical instruments.
- 5) Avoid cleaning instruments with steel wool; if want to clean with hand, use detergent soap and a soft brush and finally treat with distilled water.
- 6) You should check your autoclave; Iron, calcium, sodium or copper in water can cause rusting, staining or corrosion.
- 7) Leaving of instruments for a long time in the following solutions can also cause corrosion.
 - i) Phenol or Iodine.
 - ii) Aluminium, Barium or Calcium.
 - iii) Any Acid, Mercury or Potassium solution.
- 8) Try your utmost not to sterilize stainless steel instrument with chrome plated brass instruments, as any break in chrome surface may cause pitting or rusty surface in instruments.
- 9) Muslin wraps should be rinsed sufficiently by the laundry, so that all chlorine bleach may remove, which can cause corrosion or custing during autoclaving.

We hope that by acting on the above instructions, your problems relating to rusting, staining or pitting would abolish and the age of your instrument would increase.

MANAGEMENT

GENERAL OPERATING SCISSORS





STANDARD OPERATING SCISSORS

Sharp/Sharp Straight

HI-04-001	11.5 cm, 4 1/2"
HI-04-002	13.0 cm, 5"
HI-04-003	14.0 cm, 5 1/2"
HI-04-004	15.0 cm, 6"
HI-04-005	16.0 cm, 6 1/4"
HI-04-006	17.5 cm, 7"
HI-04-007	18.5 cm, 7 1/4"
HI-04-008	20.0 cm, 8"



STANDARD OPERATING SCISSORS

Sharp/Sharp Curved

HI-04-009	11.5 cm, 4 1/2"
HI-04-010	13.0 cm, 5"
HI-04-011	14.0 cm, 5 1/2"
HI-04-012	15.0 cm, 6"
HI-04-013	16.0 cm, 6 1/4"
HI-04-014	17.5 cm, 7"
HI-04-015	18.5 cm, 7 1/4"
HI-04-016	20.0 cm, 8"



STANDARD OPERATING SCISSORS

Blunt/Blunt Straight

HI-04-017	11.5 cm, 4 1/2"
HI-04-018	13.0 cm, 5"
HI-04-019	14.0 cm, 5 1/2"
HI-04-020	15.0 cm, 6"
HI-04-021	16.0 cm, 6 1/4"
HI-04-022	17.5 cm, 7"
HI-04-023	18.5 cm, 7 1/4"
HI-04-024	20.0 cm, 8"



STANDARD OPERATING SCISSORS

Blunt/Blunt Curved

HI-04-025	11.5 cm, 4 1/2"
HI-04-026	13.0 cm, 5"
HI-04-027	14.0 cm, 5 1/2"
HI-04-028	15.0 cm, 6"
HI-04-029	16.0 cm, 6 1/4"
HI-04-030	17.5 cm, 7"
HI-04-031	18.5 cm, 7 1/4"
HI-04-032	20.0 cm, 8"



STANDARD OPERATING SCISSORS

Sharp/Blunt Straight

HI-04-033	11.5 cm, 4 1/2"
HI-04-034	13.0 cm, 5"
HI-04-035	14.0 cm, 5 1/2"
HI-04-036	15.0 cm, 6"
HI-04-037	16.0 cm, 6 1/4"
HI-04-038	17.5 cm, 7"
HI-04-039	18.5 cm, 7 1/4"
HI-04-040	20.0 cm, 8"



STANDARD OPERATING SCISSORS

Sharp/Blunt Curved

HI-04-041	11.5 cm, 4 1/2"
HI-04-042	13.0 cm, 5"
HI-04-043	14.0 cm, 5 1/2"
HI-04-044	15.0 cm, 6"
HI-04-045	16.0 cm, 6 1/4"
HI-04-046	17.5 cm, 7"
HI-04-047	18.5 cm, 7 1/4"
HI-04-048	20.0 cm, 8"

OPERATING & DISSECTING SCISSORS



Bedienen & Präparierscheren

D'exploitation & Ciseaux de dissection

Funcionamiento y disección Tijeras

Forbici di funzionamento e di dissezione



HI-04-049
12.5 cm, 5"
B/B Str.

**LANDOLT OPERATING
SCISSORS**



HI-04-050
12.5 cm, 5"
B/B Cvd.



**DEAVER OPERATING,
DISSECTING SCISSORS**



HI-04-051
14.5 cm, 5 3/4"
S/B Str.



HI-04-052
14.5 cm, 5 3/4"
S/B Cvd.



HI-04-053
14.5 cm, 5 3/4"
S/B Str.
Saw Edge.



HI-04-054
14.5 cm, 5 3/4"
S/B Cvd.
Saw Edge.



HI-04-055
14.5 cm, 5 3/4"
S/B Str.

**LONGE OPERATING,
DISSECTING SCISSORS**



HI-04-056
14.5 cm, 5 3/4"
S/B Cvd.



HI-04-057
14.5 cm, 5 3/4"
S/B Str.
Saw Edge.



HI-04-058
14.5 cm, 5 3/4"
S/S Cvd.
Saw Edge.



INCISION OPERATING SCISSORS

1/1 original
5"

HI-04-059 13.0 cm, 5"
HI-04-060 14.5 cm, 5 3/4"
S/S Str.



MAYO DISSECTING & GYNECOLOGICAL SCISSORS



1/1 original
5 3/4"

1/4

- HI-04-061**
14.5 cm 5 3/4"
- HI-04-062**
15.0 cm 6"
- HI-04-063**
17.0 cm 6 3/4"
- HI-04-064**
19.5 cm 7 3/4"
- HI-04-065**
21.0 cm 8 1/4"
- HI-04-066**
23.0 cm 9"
B/B, Str.



1/4

- HI-04-067**
14.5 cm 5 3/4"
- HI-04-068**
15.0 cm 6"
- HI-04-069**
17.0 cm 6 3/4"
- HI-04-070**
19.5 cm 7 3/4"
- HI-04-071**
21.0 cm 8 1/4"
- HI-04-072**
23.0 cm 9"
B/B, Cvd.



MAYO-STILTE DISSECTING & GYNECOLOGICAL SCISSORS



1/1 original
6"

1/4

- HI-04-073**
15.0 cm 6"
- HI-04-074**
17.0 cm 6 3/4"
- HI-04-075**
19.5 cm 7 3/4"
- HI-04-076**
21.0 cm 8 3/4"
B/B, Str.



1/4

- HI-04-077**
15.0 cm 6"
- HI-04-078**
17.0 cm 6 3/4"
- HI-04-079**
19.5 cm 7 3/4"
- HI-04-080**
21.0 cm 8 3/4"
B/B, Cvd.



MAYO-HARRINGTON DISSECTING & GYNECOLOGICAL SCISSORS



1/4

- HI-04-081**
22.5 cm 9"
B/B, Str.



1/4

- HI-04-082**
22.5 cm 9"
B/B, Cvd.



MAYO-NOBLE DISSECTING & GYNECOLOGICAL SCISSORS



1/4

- HI-04-083**
17.0 cm, 6 3/4"
B/B, Cvd.



1/4

- HI-04-084**
17.0 cm, 6 3/4"
B/B, Str.

OPERATING & DISSECTING SCISSORS



Bedienen & Präparierscheren

D'exploitation & Ciseaux de dissection

Funcionamiento y disección Tijeras

Forbici di funzionamento e di dissezione



**METZENBAUM OPERATING,
DISSECTING SCISSORS
BLUNT/BLUNT STR.**



HI-04-085
14.5 cm 5 3/4"
B/B Str.



HI-04-086
18 cm 7"
B/B Str.



HI-04-087
20 cm 8"
B/B Str.



HI-04-088
23 cm 9"
B/B Str.



HI-04-089
25 cm 10"
B/B Str.



HI-04-090
28 cm 11"
B/B Str.



**METZENBAUM OPERATING,
DISSECTING SCISSORS
BLUNT/BLUNT CVD.**



HI-04-091
14.5 cm 5 3/4"
B/B Cvd.



HI-04-092
18 cm 7"
B/B Cvd.



HI-04-093
20 cm 8"
B/B Cvd.



HI-04-094
23 cm 9"
B/B Cvd.



HI-04-095
25 cm 10"
B/B Cvd.



HI-04-096
28 cm 11"
B/B Cvd.



**METZENBAUM-NELSON
SCISSORS**



1/1 original
7"

- HI-04-097 18.0 cm 7"
- HI-04-098 20.0 cm 8"
- HI-04-099 23.0 cm 9"
- HI-04-100 25.0 cm 9 3/4"
- HI-04-101 28.0 cm 11"
- HI-04-102 30.0 cm 12"
B/B Str.



1/1 original
7"

- HI-04-103 18.0 cm 7"
- HI-04-104 20.0 cm 8"
- HI-04-105 23.0 cm 9"
- HI-04-106 25.0 cm 9 3/4"
- HI-04-107 28.0 cm 11"
- HI-04-108 30.0 cm 12"
S/B Str.



- HI-04-109 18.0 cm 7"
- HI-04-110 20.0 cm 8"
- HI-04-111 23.0 cm 9"
- HI-04-112 25.0 cm 9 3/4"
- HI-04-113 28.0 cm 11"
- HI-04-114 30.0 cm 12"
S/S Str.



**METZENBAUM-NELSON
SCISSORS**



1/1 original
7"

- HI-04-115 18.0 cm 7"
- HI-04-116 20.0 cm 8"
- HI-04-117 23.0 cm 9"
- HI-04-118 25.0 cm 9 3/4"
- HI-04-119 28.0 cm 11"
- HI-04-120 30.0 cm 12"
B/B Cvd.



1/1 original
7"

- HI-04-121 18.0 cm 7"
- HI-04-122 20.0 cm 8"
- HI-04-123 23.0 cm 9"
- HI-04-124 25.0 cm 9 3/4"
- HI-04-125 28.0 cm 11"
- HI-04-126 30.0 cm 12"
S/B Cvd.



- HI-04-127 18.0 cm 7"
- HI-04-128 20.0 cm 8"
- HI-04-129 23.0 cm 9"
- HI-04-130 25.0 cm 9 3/4"
- HI-04-131 28.0 cm 11"
- HI-04-132 30.0 cm 12"
S/S Cvd.

OPERATING & DISSECTING SCISSORS



Bedienen & Präparierscheren

D'exploitation & Ciseaux de dissection

Funcionamiento y disección Tijeras

Forbici di funzionamento e di dissezione



**METZENBAUM-FINO
SCISSORS**



HI-04-133
14.5 cm 5 3/4"
B/B Str.



HI-04-134
14.5 cm 5 3/4"
S/S Str.



HI-04-135
14.5 cm 5 3/4"
B/B Cvd.



HI-04-136
14.5 cm 5 3/4"
S/S Cvd.



**METZENBAUM-FINO
SCISSORS**



HI-04-137
18 cm 7"
HI-04-138
220 cm 8"
HI-04-139
230 cm 9"
B/B Str.

1/1 original
7"



HI-04-140
18 cm 7"
HI-04-141
220 cm 8"
HI-04-142
230 cm 9"
B/B Cvd.

1/1 original
7"



HI-04-143
18 cm 7"
HI-04-144
220 cm 8"
HI-04-145
230 cm 9"
B/B Slightly
Cvd.



**RICHTER OPERATING,
DISSECTING SCISSORS**


HI-04-146 14.5 cm, 5 3/4"
Angled



**JOSEPH OPERATING,
DISSECTING SCISSORS**



HI-04-147
14.0 cm, 5 1/2"
S/S Str.



HI-04-148
14.0 cm, 5 1/2"
S/S Cvd



**KELLY OPERATING,
DISSECTING SCISSORS**



HI-04-149
16.0 cm, 6 1/4"
HI-04-150
18.0 cm, 7"
S/S Str.

1/1 original
6 1/4"



HI-04-151
16.0 cm, 6 1/4"
HI-04-152
18.0 cm, 7"
S/S Cvd.



**DOYEN OPERATING,
DISSECTING SCISSORS**

HI-04-153 17.5 cm, 7"
B/B, Str.



**DOYEN OPERATING,
DISSECTING SCISSORS**


HI-04-154 17.5 cm, 7"
B/B, Cvd.




 1/1 original
 8"

 1/1 original
 8"

 1/1 original
 8"

 1/1 original
 8"


HI-04-157 20.0 cm, 8"	HI-04-159 20.0 cm, 8"	HI-04-161 20.0 cm, 8"	HI-04-163 20.0 cm, 8"	HI-04-165 20.0 cm, 8"
HI-04-158 23.0 cm, 9"	HI-04-160 23.0 cm, 9"	HI-04-162 23.0 cm, 9"	HI-04-164 23.0 cm, 9"	HI-04-166 23.0 cm, 9"
Str, S/B	Str, S/S	Cvd, B/B	Cvd, S/B	Cvd, S/S

SIMS DISSECTING & GYNECOLOGICAL SCISSORS
HI-04-155 20.0 cm, 8"
HI-04-156 23.0 cm, 9"
 Str, B/B

Iris Scissors

 1/1 original
 3 1/2"


HI-04-167 09 cm, 3 1/2"	HI-04-171 09 cm, 3 1/2"
HI-04-168 10.5 cm, 4 1/4"	HI-04-172 10.5 cm, 4 1/4"
HI-04-169 11.5 cm, 4 1/2"	HI-04-173 11.5 cm, 4 1/2"
HI-04-170 12.0 cm, 4 3/4"	HI-04-174 12.0 cm, 4 3/4"
S/S Str.	S/S Cvd.



HI-04-175
 11.5 cm, 4 1/2"
 Angular Left

HI-04-176
 11.5 cm, 4 1/2"
 Laterally Cvd.

HI-04-177
 11.5 cm, 4 1/2"
 Angular Right



VASCULAR SCISSORS
HI-04-178 19.0 cm, 7 1/2"
S/S Angled 125°



COOLEY REVERSE CUT SCISSORS
HI-04-179 17.1 cm, 6 3/4"
S/S angled 125°



HI-04-180
19.0 cm, 7 1/2"
25°



HI-04-181
19.0 cm, 7 1/2"
45°



HI-04-182
19.0 cm, 7 1/2"
60°



HI-04-183
19.0 cm, 7 1/2"
125°



HEGEMANN-DIETRICH VASCULAR SCISSORS
S/S Angled




25°
HI-04-184 16.0 cm, 6 1/4"
HI-04-185 19.0 cm, 7 1/2"



45°
HI-04-186 16.0 cm, 6 1/4"
HI-04-187 19.0 cm, 7 1/2"



60°
HI-04-188 16.0 cm, 6 1/4"
HI-04-189 19.0 cm, 7 1/2"



90°
HI-04-190 16.0 cm, 6 1/4"
HI-04-191 19.0 cm, 7 1/2"



125°
HI-04-192 16.0 cm, 6 1/4"
HI-04-193 19.0 cm, 7 1/2"

VASCULAR SCISSORS



Gefäßschere | Ciseaux vasculaires
Tijeras vasculares | Forbici vascolari



DE BAKEY VASCULAR SCISSORS
S/S Angled




1/1 original
6 1/4"

25°
HI-04-194 16.0 cm, 6 1/4"
HI-04-195 19.0 cm, 7 1/2"
HI-04-196 23.0 cm, 9"
HI-04-197 28.0 cm, 11"



1/1 original
6 1/4"

45°
HI-04-198 16.0 cm, 6 1/4"
HI-04-199 19.0 cm, 7 1/2"
HI-04-200 23.0 cm, 9"
HI-04-201 28.0 cm, 11"




60°
HI-04-202 16.0 cm, 6 1/4"
HI-04-203 19.0 cm, 7 1/2"
HI-04-204 23.0 cm, 9"
HI-04-205 28.0 cm, 11"



COOLEY VASCULAR SCISSORS
HI-04-206 18.4 cm, 7 1/2" 1 mm
HI-04-207 18.4 cm, 7 1/2" 1.5 mm
probe point on tip.
S/S Angled 45°



O'BRIEN VASCULAR SCISSORS
HI-04-208 9.0 cm, 3 1/2"
S/S Angled




HEATH VASCULAR SCISSORS
HI-04-209 15.0 cm, 6"
S/S Double Cvd.



HI-04-210
17.0 cm, 6 3/4"
Serrated Blade S/S
Double Cvd.

**DEAN-DELLICATE
DISSECTING SCISSORS**




Locklin Dissecting Scissors
HI-04-211 16.5 cm, 6 1/2"
S/S



FAVALORO CORONARY
HI-04-212 10.5 cm, 4"
S/S Angled 90°



ENT-DEAN SCISSORS
HI-04-213 17 cm, 6 3/4"
S/S Double Cvd.
Serrated



HI-04-214
15.5 cm, 6"
S/S

HI-04-215
15.5 cm, 6"
S/S

**DE BAKKEY VASCULAR
SCISSORS**

VASCULAR SCISSORS




Gefäßschere

Ciseaux vasculaires

Tijeras vasculares

Forbici vascolari



¼

HI-04-216

10 cm 4"

HI-04-217

11.5 cm 4 1/2"

S/S Str.

STEVENS SCISSORS



¼


HI-04-218

10 cm 4"

HI-04-219

11.5 cm 4 1/2"

S/S Cvd.



¼

HI-04-220

10 cm 4"

HI-04-221

11.5 cm 4 1/2"

B/B Str.

STEVENS SCISSORS



¼

HI-04-222

10 cm 4"

HI-04-223

11.5 cm 4 1/2"

B/B Cvd.



¼

NAIL SCISSORS

HI-04-224 9 cm, 3 1/2"

S/S, Cvd.



¼

NAIL SCISSORS

HI-04-225 9 cm, 3 1/2"

S/S, Cvd.



¼

NAIL SCISSORS

HI-04-226 9 cm, 3 1/2"


Arrow Point Curved

STITCH SCISSORS


stitch Schere | point de ciseaux
 Coser Tijeras | Stitch Forbici



STITCHES SCISSORS
 HI-04-227 12.0 cm, 4 3/4"
 Str.



NORTHBENT STITCHES SCISSORS
 HI-04-228 12.5 cm, 5"
 Cvd.



LITTAUER STITCH SCISSORS
 HI-04-229 14 cm, 5 1/2"
 Str.




BUCK STITCH SCISSORS

HI-04-230
 14.5 cm, 5 3/4"
 Str.




SPENCER STITCH SCISSORS

HI-04-231
 18.0 cm, 7"
 Str.



HI-04-232
 9 cm, 3 1/2"
 Str.

HI-04-233
 11.5 cm, 4 1/2"
 Str.




SHORTBENT STITCH SCISSORS
 HI-04-234 9 cm, 3 1/2"
 Cvd.

TONSIL, NASAL SCISSORS



Tonsillensteinen Schere, Nasenschere
Tijeras de amígdalas, tijeras nasales

Ciseaux d'amygdale, Ciseaux nasaux
Tonsille forbici, forbici nasali



COTTLE SCISSORS
HI-04-235 16 cm 6 1/4"
 B/B, Angled




FOMON SCISSORS
HI-04-236 14.5 cm 5 1/2"
 B/B, Angled



SAILER SCISSORS
HI-04-237 17.0 cm 6 3/4"
 B/B, Angled



HEYMANN SCISSORS
HI-04-238 18 cm, 7"
HI-04-239 19.5 cm, 7 3/4"
 B/B, Angled



HEYMANN SCISSORS
HI-04-240 18.0 cm, 7"
 B/B, Angled



OLIVECRONA VASCULAR SCISSORS
HI-04-241 23.0 cm, 9"
 B/B

GUM & OPERATING SCISSORS



Chirurgische scheren, Gum Scheren
Goma Tijeras y tijeras de funcionamiento

Gum Ciseaux et Ciseaux d'exploitation
Gum Forbici & Forbici Operativi



¼

GOOD SCISSORS
HI-04-242 19.0 cm, 7 ½"
Probe End, Cvd.



¼

PRINCE SCISSORS
HI-04-243 17.5 cm, 7"
B/B, Cvd.



¼


HALLE SCISSORS
HI-04-244 19.0 cm, 7 ½"
B/B, Cvd.



HI-04-245
14.5 cm, 5 ¾"
Pointed Tip
S/S, Str.



HI-04-246
14.5 cm, 5 ¾"
Pointed Tip
S/S, Cvd.



¼

MOD. BONN SCISSORS
HI-04-247 8.0 cm, 3"
S/S, Str.




¼

MOD. BONN SCISSORS
HI-04-248 8.0 cm, 3"
S/S, Cvd.

WIRE CUTTING SCISSORS

Draht-Ausschnitt-Scheren
Tijeras de corte de alambre

Ciseaux de coupe de câbles
Filo da taglio Forbici



HI-04-249
10.5 cm 4 1/4"
HI-04-250
12.0 cm 4 3/4"
B/B Str.
HI-04-251
10.5 cm 4 1/4"
HI-04-252
12.0 cm 4 3/4"
B/B Str.
Saw Edge



HI-04-253
10.5 cm 4 1/4"
HI-04-254
12.0 cm 4 3/4"
S/S Str.
HI-04-255
10.5 cm 4 1/4"
HI-04-256
12.0 cm 4 3/4"
S/S Str.
Saw Edge



HI-04-257
10.0 cm 4"
HI-04-258
12.0 cm 4 3/4"
B/B Cvd.
HI-04-259
10.0 cm 4"
HI-04-260
12.0 cm 4 3/4"
B/B Cvd.
Saw Edge



HI-04-261
10.0 cm 4"
HI-04-262
12.0 cm 4 3/4"
S/S Cvd.
HI-04-263
10.0 cm 4"
HI-04-264
12.0 cm 4 3/4"
S/S Cvd.
Saw Edge

BEEBEE WIRE CUTTING SCISSORS

BEEBEE WIRE CUTTING SCISSORS



SYSTRUNK
HI-04-265
13.0 cm 5"
S/S, Str.



SYSTRUNK
HI-04-266
13.0 cm 5"
S/S, Cvd.

SYSTRUNK SCISSORS



UNIVERSAL WIRE CUTTING SCISSORS
HI-04-267 12 cm, 4 3/4"

NAIL NIPPERS

Nagelzangen | pinces à ongles
Nail Nippers | Nail Tronchesi



NAIL NIPPERS
HI-04-268 11 cm 4 1/4"
HI-04-269 13 cm 5"
HI-04-270 15 cm 6"



NAIL NIPPERS
HI-04-271 14 cm 5 1/2"



NAIL NIPPERS
HI-04-272 14 cm 5 1/2"



TURNBULL NIPPERS
HI-04-273 11.5 cm 4 1/2"
HI-04-274 13 cm 5"



NAIL NIPPERS
HI-04-275 10.5 cm 4 1/4"



COHEN NIPPERS
HI-04-276 10.5 cm 4 1/4"

MICRO SCISSORS

Mikroscheren | micro Ciseaux
 Micro Tijeras | Micro Forbici




ARTERY SCISSORS
 Probe end Angled




HI-04-277
 18.0 cm, 7"
 25°




HI-04-278
 18.0 cm, 7"
 45°



HI-04-279
 18.0 cm, 7"
 60°



HI-04-280
 18.0 cm, 7"
 90°



HI-04-281
 18.0 cm, 7"
 125°



HI-04-282
 18.0 cm, 7"
 25°



HI-04-283
 18.0 cm, 7"
 45°



HI-04-284
 18.0 cm, 7"
 60°



HI-04-285
 18.0 cm, 7"
 90°



HI-04-286
 18.0 cm, 7"
 125°

MICRO SCISSORS

Mikroscheren | *micro Ciseaux*
 Micro Tijeras | *Micro Forbici*



¼

HI-04-287
 9.0 cm, 3 ½"
 S/B Str.



¼

HI-04-288
 9.0 cm, 3 ½"
 S/B Cvd.



¼

HI-04-289
 9.0 cm, 3 ½"
 S/S Str.



¼

HI-04-290
 9.0 cm, 3 ½"
 S/S Cvd.



¼

HI-04-291
 11.4 cm, 4 ½"
 B/B Str.



¼

HI-04-292
 11.4 cm, 4 ½"
 B/B Cvd.



¼

HI-04-293
 11.4 cm, 4 ½"
 S/B Str.



¼

HI-04-294
 11.4 cm, 4 ½"
 S/B Cvd.



¼

HI-04-295
 11.4 cm, 4 ½"
 S/S Str.




¼

HI-04-296
 11.4 cm, 4 ½"
 S/S Cvd.



¼

HI-04-297
 14 cm, 5 ½"
 B/B Str.



¼

HI-04-298
 14 cm, 5 ½"
 B/B Cvd.



¼

HI-04-299
 14 cm, 5 ½"
 S/B Str.



¼

HI-04-300
 14 cm, 5 ½"
 S/B Cvd.



¼

HI-04-301
 14 cm, 5 ½"
 S/S Str.



¼

HI-04-302
 14 cm, 5 ½"
 S/S Cvd.



CASTROVIEJO SCISSORS
 Spring Type

MICRO SCISSORS

Mikroscheren | micro Ciseaux
 Micro Tijeras | Micro Forbici



NOYES SCISSORS
 Spring Handle
HI-04-306 12 cm, 4 3/4"
 Angled,



WECKER SCISSORS
 Angled



BARRAQUER SCISSORS
HI-04-310 Sharp Points
HI-04-311 Blunt Points
 5.5 cm, 2 1/4"



CASTROVIEJO CORNEAL MICRO SCISSORS
HI-04-312 11 cm, 4 1/3"
 Left Angled,
 Spring Handle



CASTROVIEJO CORNEAL MICRO SCISSORS
HI-04-313 11 cm, 4 1/3"
 Right Angled,
 Spring Handle





TROUTMAN MICRO SCISSORS
HI-04-314 10.0 cm, 4"
 Left Angled,
 Spring Handle



TROUTMAN MICRO SCISSORS
HI-04-315 10.0 cm, 4"
 Right Angled,
 Spring Handle



WESTCOTT MICRO SCISSORS



MICRO SCISSORS

Mikroscheren | *micro Ciseaux*
 Micro Tijeras | *Micro Forbici*



MICRO SCISSORS
HI-04-320 12.0 cm, 4 3/4 S/S, Angled
 S/S, Cvd.



MICRO SCISSORS



1/4
HI-04-321
 18.0 cm, 7"
 25°



1/4
HI-04-322
 18.0 cm, 7"
 45°



1/4
HI-04-323
 18.0 cm, 7"
 60°



1/4
HI-04-324
 18.0 cm, 7"
 90°



1/4
HI-04-325
 18.0 cm, 7"
 125°



MICRO SCISSORS
 S/S, Angled pointed probe



1/4
HI-04-326
 18.0 cm, 7"
 25°



1/4
HI-04-327
 18.0 cm, 7"
 45°



1/4
HI-04-328
 18.0 cm, 7"
 60°



1/4
HI-04-329
 18.0 cm, 7"
 90°



1/4
HI-04-330
 18.0 cm, 7"
 125°



YASARGILL SCISSORS



1/4
HI-04-331 16 cm, 6 1/4"
HI-04-332 20 cm, 8"
HI-04-333 22.5 cm, 9"
 S/S, Str.



1/4 P
HI-04-334 16 cm, 6 1/4"
HI-04-335 20 cm, 8"
HI-04-336 22 cm, 8 3/4"
 S/S, Cvd.



1/4
HI-04-337 16 cm, 6 1/4"
HI-04-338 20 cm, 8"
HI-04-339 22 cm, 8 3/4"
 B/B, Str.



1/4
HI-04-340 16 cm, 6 1/4"
HI-04-341 20 cm, 8"
HI-04-342 22 cm, 8 3/4"
 S/S, Str.



1/4
HI-04-343 16 cm, 6 1/4"
HI-04-344 20 cm, 8"
HI-04-345 22 cm, 8 3/4"
 S/S, Cvd.



1/4
HI-04-346 16 cm, 6 1/4"
HI-04-347 20 cm, 8"
HI-04-348 22 cm, 8 3/4"
 B/B, Cvd.



NOYES MICRO SCISSORS
HI-04-349 15.0 cm, 6"
 S/S



NOYES MICRO SCISSORS
HI-04-350 14.5 cm, 5 3/4"
 S/S

UNIVERSAL SCISSORS

Universal-Schere | Ciseaux universelles
Tijeras universales | Forbici universali



UNIVERSAL PARAMEDIC SHEAR

HI-04-351 14 cm 5 3/4" (BLACK)
HI-04-352 14 cm 5 3/4" (BLUE)
HI-04-353 14 cm 5 3/4" (ROYAL BLUE)
HI-04-354 14 cm 5 3/4" (BURGANDY)
HI-04-355 14 cm 5 3/4" (GREEN)
HI-04-356 14 cm 5 3/4" (MAGENTA)
HI-04-357 14 cm 5 3/4" (ORANGE)
HI-04-358 14 cm 5 3/4" (PINK)
HI-04-359 14 cm 5 3/4" (PURPLE)
HI-04-360 14 cm 5 3/4" (RED)
HI-04-361 14 cm 5 3/4" (TEAL)
HI-04-362 14 cm 5 3/4" (WHITE)
HI-04-363 14 cm 5 3/4" (YELLOW)

HI-04-364 16 cm 6 1/4" (BLACK)
HI-04-365 16 cm 6 1/4" (BLUE)
HI-04-366 14 cm 5.5" (ROYAL BLUE)
HI-04-367 14 cm 5.5" (BURGANDY)
HI-04-368 16 cm 6 1/4" (GREEN)
HI-04-369 16 cm 6 1/4" (MAGENTA)
HI-04-370 16 cm 6 1/4" (ORANGE)
HI-04-371 16 cm 6 1/4" (PINK)
HI-04-372 16 cm 6 1/4" (PURPLE)
HI-04-373 16 cm 6 1/4" (RED)
HI-04-374 16 cm 6 1/4" (TEAL)
HI-04-375 16 cm 6 1/4" (WHITE)
HI-04-376 16 cm 6 1/4" (YELLOW)

143 C/290 F Autoclavable

**UNIVERSAL
PARAMEDIC SHEAR
HI-04-377**
14 cm 5.5"
(GREEN)



**UNIVERSAL
PARAMEDIC SHEAR
HI-04-382**
14 cm 5.5"
(PINK)



**UNIVERSAL
PARAMEDIC SHEAR
HI-04-381**
14 cm 5.5"
(BURGANDY)



**UNIVERSAL
PARAMEDIC SHEAR
HI-04-378**
14 cm 5.5"
(ROYAL BLUE)



**UNIVERSAL
PARAMEDIC SHEAR
HI-04-379**
14 cm 5.5"
(RED)



**UNIVERSAL
PARAMEDIC SHEAR
HI-04-380**
14 cm 5.5"
(BLACK)

UNIVERSAL SCISSORS

Universal-Schere | Ciseaux universelles
Tijeras universales | Forbici universali



**UNIVERSAL PARAMEDIC
SHEAR PROB END
HI-04-383
14 cm 5.5"
(YELLOW)**



**UNIVERSAL PARAMEDIC
SHEAR PROB END
HI-04-390
14 cm 5.5"
(ROYAL BLUE)**



**UNIVERSAL PARAMEDIC
SHEAR PROB END
HI-04-384
14 cm 5.5"
(Chartreuse)**



**UNIVERSAL PARAMEDIC
SHEAR PROB END
HI-04-389
14 cm 5.5"
(WHITE)**



**UNIVERSAL PARAMEDIC
SHEAR PROB END
HI-04-385
14 cm 5.5"
(RED)**



**UNIVERSAL PARAMEDIC
SHEAR PROB END
HI-04-386
14 cm 5.5"
(BLACK)**



**UNIVERSAL PARAMEDIC
SHEAR PROB END
HI-04-388
14 cm 5.5"
(GREEN)**

**UNIVERSAL PARAMEDIC
SHEAR PROB END
HI-04-387
14 cm 5.5"
(BURGUNDY)**

UNIVERSAL SCISSORS

Universal-Schere | Ciseaux universelles
Tijeras universales | Forbici universali

**UNIVERSAL
PARAMEDIC SHEAR**
HI-04-391
14 cm 5.5"
(BLACK)



**UNIVERSAL
PARAMEDIC SHEAR**
HI-04-395
14 cm 5.5"
(BLUE)



**UNIVERSAL
PARAMEDIC SHEAR**
HI-04-392
14 cm 5.5"
(ORANGE)



**UNIVERSAL
PARAMEDIC SHEAR**
HI-04-396
14 cm 5.5"
(YELLOW)



**UNIVERSAL
PARAMEDIC SHEAR**
HI-04-393
14 cm 5.5"
(WHITE)



**UNIVERSAL
PARAMEDIC SHEAR**
HI-04-397
14 cm 5.5"
(BURGUNDY)



**UNIVERSAL
PARAMEDIC SHEAR**
HI-04-394
14 cm 5.5"
(ORANGE)



**UNIVERSAL
PARAMEDIC SHEAR**
HI-04-398
14 cm 5.5"
(DEEP YELLOW)



PLASTER SCISSORS



Gips Schere | plâtre Ciseaux
Tijeras Yeso | intonaco Forbici



ESMARCH KNIFE
HI-04-399 18 cm, 7"



BRUNS SCISSORS



HI-04-400
24 cm, 9 1/2"
Cvd, Probe End



HI-04-401
24 cm, 9 1/2"
Cvd, Serrated
Probe End



**LISTER-EXCENTRIC
PLASTER SCISSORS**
HI-04-402 16 cm, 6 1/4"
HI-04-403 19 cm, 7 1/2"
Angled, Probe End



Edge Upper
is Serrated



**LISTER BANDAGE PLASTER
Scissors, one large ring**
HI-04-404 9 cm, 3 1/2"
HI-04-405 11 cm, 4 1/3"
HI-04-406 14 cm, 5 1/2"
HI-04-407 15 cm, 6"
HI-04-408 18 cm, 7"
HI-04-409 20 cm, 8"
Angled, Probe End



BERGMANN SCISSORS
HI-04-410 23 cm, 9"
Angled, Probe End



ESMARCH SCISSORS
HI-04-411 20 cm, 8"
HI-04-412 23 cm, 9"
Angled, Probe End



**LISTER BANDAGE
PLASTER SCISSORS**
HI-04-413 9 cm, 3 1/2"
HI-04-414 11 cm, 4 1/3"
HI-04-415 14 cm, 5 1/2"
HI-04-416 15 cm, 6"
HI-04-417 18 cm, 7"
HI-04-418 20 cm, 8"
Angled, Probe End



PLASTER SCISSORS
HI-04-419 17 cm, 6 3/4"
Angled



**BURNHAM BANDAGE
PLASTER SCISSORS**
HI-04-420 14 cm, 5 1/2"
Angled

PLASTER SCISSORS

Gips Schere | plâtre Ciseaux
Tijeras Yeso | intonaco Forbici



BANDAGE SCISSORS
HI-04-421
16.0 cm 6 1/4"



**SMITH (Mod. Usa)
PLASTER SCISSORS**
HI-04-422 16 cm, 6 1/4"
HI-04-423 18 cm, 7"
HI-04-424 20 cm, 8"
HI-04-425 23 cm, 9"



**LORENZ PLASTER
SCISSORS**
HI-04-426 23 cm, 9"
B/B, Cvd



**LANGE BANDAGE
SCISSORS**
HI-04-427 18.0 cm 7"
B/B



**BRAUN-STADLER
PLASTER SCISSORS**
HI-04-428 14 cm, 5 1/2"
HI-04-429 22 cm, 8 3/4"
B/B, Angled



**BUSCH PLASTER
SCISSORS**
HI-04-430 13 cm, 5"
HI-04-431 16 cm, 6 1/4"
B/B, Cvd.



**KNOWLES PLASTER
SCISSORS**
HI-04-432 14 cm, 5 1/2"
Str. Probe End



**KNOWLES PLASTER
SCISSORS**
HI-04-433 14 cm, 5 1/2"
Angled, Probe End



**WALDMANN PLASTER
SCISSORS**
HI-04-434 18 cm, 7"
Probe End

PLASTER SCISSORS

Gips Schere | plâtre Ciseaux
Tijeras Yeso | intonaco Forbici



WALDMANN SCISSORS
HI-04-435 18 cm, 7"



SCHUMACHER PLASTER SCISSORS
HI-04-436 16 cm, 6 1/4"



MODIFIED-USA PLASTER SCISSORS
HI-04-437 10.5 cm, 4 1/4"



MARBACH PLASTER SCISSORS
HI-04-438 15.5 cm, 6"
Probe End, Serrated blade



SEUTIN PLASTER SCISSORS
HI-04-439 20.0 cm 8"
Probe End



SEUTIN PLASTER SCISSORS
HI-04-440 23.0 cm 9"
Probe End

PLASTER SHEARS

GipsScheren | plâtres Ciseaux
Tijeras Yeso | Cesoie in gesso



STILLE PLASTER SHEARS
HI-04-441 23 cm, 9"
HI-04-442 Spare Blade



STILLE PLASTER SHEARS
HI-04-443 26 cm, 10 1/2"
HI-04-444 Spare Blade



STILLE PLASTER SHEARS
HI-04-445 37 cm, 14 1/2"
HI-04-446 Spare Blade



STILLE PLASTER SHEARS
HI-04-447 24 cm, 9 1/2"
HI-04-448 Spare Blade

PLASTER SHEARS


GipsScheren | plâtres Ciseaux
Tijeras Yeso | Cesoie in gesso



HENNIG PLASTER CAST SPREADER
HI-04-449 28 cm, 11"



MODEL USA CAST SPREADER PLASTER
HI-04-450 23 cm, 9" 3 Prong



SPREADER PLASTER
HI-04-451 29 cm, 11 1/2"



WOLFF SPREADER PLASTER
HI-04-452 18 cm, 7"
HI-04-453 24 cm, 9 1/2"



ENGEL PLASTER SAW
HI-04-454 14 cm, 5 1/2"



CHARRIERE PLASTER SAW
HI-04-455 25 cm, 10"
HI-04-456 27 cm, 10 3/4"
Double saw

SUPER CUT SCISSORS

Super Cut Schere | Ciseaux SuperCut
 Súper Cut Tijeras | Forbici SuperCut



STANDARD OPERATING
 SUPER CUT SCISSORS




HI-04-457 14.5 cm 5 3/4"
HI-04-458 16.5 cm 6 1/2"
 Str. B/B

HI-04-459 14.5 cm 5 3/4"
HI-04-460 16.5 cm 6 1/2"
 Str. S/B

HI-04-461 14.5 cm 5 3/4"
HI-04-462 16.5 cm 6 1/2"
 Str. S/S



STANDARD OPERATING
 SUPER CUT SCISSORS



HI-04-463 14.5 cm 5 3/4"
HI-04-464 16.5 cm 6 1/2"
 Cvd. B/B

HI-04-465 14.5 cm 5 3/4"
HI-04-466 16.5 cm 6 1/2"
 Cvd. S/B

HI-04-467 14.5 cm 5 3/4"
HI-04-468 16.5 cm 6 1/2"
 Cvd. S/S



LAXER OPERATING SUPER
 CUT SCISSORS



HI-04-469 16 cm 6 1/4"
HI-04-470 21 cm 8 1/4"
 B/B, Str.



HI-04-471 16.5 cm 5 1/2"
HI-04-472 21 cm 8 1/4"
 B/B, Cvd.

SUPER CUT SCISSORS

Super Cut Schere | Ciseaux SuperCut
Súper Cut Tijeras | Forbici SuperCut



SAFETY SUPER CUT SCISSORS
HI-04-473 Str.
HI-04-474 Cvd.
9.0 cm, 3 1/2", Probe Ends



**GOLDMAN-FOX
SUPER CUT SCISSORS**
HI-04-475 Str.
HI-04-476 Cvd.
9.0 cm, 3 1/2", S/S
13 cm 5"



**STEVENS TENOTOMY
SUPER CUT SCISSORS**
HI-04-477 Str.
HI-04-478 Cvd.
9.0 cm, 3 1/2", B/B



**STEVENS TENOTOMY
SUPER CUT SCISSORS**
HI-04-479 Str.
HI-04-480 Cvd.
11.5 cm, 4 1/2", S/S



**BABY METZENBAUM
SUPER CUT SCISSORS**
HI-04-481 11.5 cm, 4 1/2"
B/B, Str.



**BABY METZENBAUM
SUPER CUT SCISSORS**
HI-04-482 11.5 cm, 4 1/2"
B/B, Cvd.



IRIS SUPER CUT SCISSORS
HI-04-483 9.0 cm, 3 1/2"
HI-04-484 11.5 cm, 4 1/2"
S/S, Str.



IRIS SUPER CUT SCISSORS
HI-04-485 9.0 cm, 3 1/2"
HI-04-486 11.5 cm, 4 1/2"
S/S, Cvd.


**Operating Scissors, S/S Str.
T/C Insert**

HI-04-487	11.5 cm, 4 1/2"
HI-04-488	13.0 cm, 5"
HI-04-489	14.0 cm, 5 1/2"
HI-04-490	15.0 cm, 6"
HI-04-491	16.0 cm, 6 1/4"
HI-04-492	17.5 cm, 7"
HI-04-493	18.5 cm, 7 1/4"
HI-04-494	20.0 cm, 8"


**Operating Scissors, S/S Cvd.
T/C Insert**

HI-04-495	11.5 cm, 4 1/2"
HI-04-496	13.0 cm, 5"
HI-04-497	14.0 cm, 5 1/2"
HI-04-498	15.0 cm, 6"
HI-04-499	16.0 cm, 6 1/4"
HI-04-500	17.5 cm, 7"
HI-04-501	18.5 cm, 7 1/4"
HI-04-502	20.0 cm, 8"


**Operating Scissors, B/B Str.
T/C Insert**

HI-04-503	11.5 cm, 4 1/2"
HI-04-504	13.0 cm, 5"
HI-04-505	14.0 cm, 5 1/2"
HI-04-506	15.0 cm, 6"
HI-04-507	16.0 cm, 6 1/4"
HI-04-508	17.5 cm, 7"
HI-04-509	18.5 cm, 7 1/4"
HI-04-510	20.0 cm, 8"


**Operating Scissors, B/B Cvd.
T/C Insert**

HI-04-511	11.5 cm, 4 1/2"
HI-04-512	13.0 cm, 5"
HI-04-513	14.0 cm, 5 1/2"
HI-04-514	15.0 cm, 6"
HI-04-515	16.0 cm, 6 1/4"
HI-04-516	17.5 cm, 7"
HI-04-517	18.5 cm, 7 1/4"
HI-04-518	20.0 cm, 8"


**Operating Scissors, S/B
Standard Str. T/C Insert**

HI-04-519	11.5 cm, 4 1/2"
HI-04-520	13.0 cm, 5"
HI-04-521	14.0 cm, 5 1/2"
HI-04-522	15.0 cm, 6"
HI-04-523	16.0 cm, 6 1/4"
HI-04-524	17.5 cm, 7"
HI-04-525	18.5 cm, 7 1/4"
HI-04-526	20.0 cm, 8"


**Operating Scissors, S/B
Standard Cvd. T/C Insert**

HI-04-527	11.5 cm, 4 1/2"
HI-04-528	13.0 cm, 5"
HI-04-529	14.0 cm, 5 1/2"
HI-04-530	15.0 cm, 6"
HI-04-531	16.0 cm, 6 1/4"
HI-04-532	17.5 cm, 7"
HI-04-533	18.5 cm, 7 1/4"
HI-04-534	20.0 cm, 8"

OPERATING SCISSORS T.C.



Chirurgische scheren
Tijeras quirúrgicas

Ciseaux chirurgicaux
Forbici per chirurgia



Metzenbaum Operating Scissors T/C Insert



HI-04-535
11.5 cm 4 1/2"
HI-04-536
14.5 cm 5 3/4"
B/B Str.



HI-04-537
11.5 cm 4 1/2"
HI-04-538
14.5 cm 5 3/4"
S/B Str.



HI-04-539
11.5 cm 4 1/2"
HI-04-540
14.5 cm 5 3/4"
S/S Str.



Metzenbaum Operating Scissors T/C Insert



HI-04-541
11.5 cm 4 1/2"
HI-04-542
14.5 cm 5 3/4"
B/B Cvd.



HI-04-543
11.5 cm 4 1/2"
HI-04-544
14.5 cm 5 3/4"
S/B Cvd.



HI-04-545
11.5 cm 4 1/2"
HI-04-546
14.5 cm 5 3/4"
S/S Cvd.



Metzenbaum-Nelson Operating Scissors T/C Insert



HI-04-547
18 cm 7"
HI-04-548
20 cm 8"
HI-04-549
23 cm 9"
HI-04-550
25 cm 10"
HI-04-551
28 cm 11"
HI-04-552
30 cm 12"
B/B Str.



HI-04-553
18 cm 7"
HI-04-554
20 cm 8"
HI-04-555
23 cm 9"
HI-04-556
25 cm 10"
HI-04-557
28 cm 11"
HI-04-558
30 cm 12"
S/B Str.



HI-04-559
18 cm 7"
HI-04-560
20 cm 8"
HI-04-561
23 cm 9"
HI-04-562
25 cm 10"
HI-04-563
28 cm 11"
HI-04-564
30 cm 12"
S/S Str.



Metzenbaum-Nelson Operating Scissors T/C Insert



HI-04-565
18 cm 7"
HI-04-566
20 cm 8"
HI-04-567
23 cm 9"
HI-04-568
25 cm 10"
HI-04-569
28 cm 11"
HI-04-570
30 cm 12"
B/B Cvd.



HI-04-571
18 cm 7"
HI-04-572
20 cm 8"
HI-04-573
23 cm 9"
HI-04-574
25 cm 10"
HI-04-575
28 cm 11"
HI-04-576
30 cm 12"
S/B Cvd.



HI-04-577
18 cm 7"
HI-04-578
20 cm 8"
HI-04-579
23 cm 9"
HI-04-580
25 cm 10"
HI-04-581
28 cm 11"
HI-04-582
30 cm 12"
S/S Cvd.



$\frac{1}{4}$
HI-04-583 14.5 cm 5 3/4"
HI-04-584 18 cm 7"
HI-04-585 20 cm 8"
HI-04-586 23 cm 9"
 B/B Str.

Metzenbaum-Fino
Operating Scissors
T/C Insert



$\frac{1}{4}$
HI-04-587 14.5 cm 5 3/4"
HI-04-588 18 cm 7"
HI-04-589 20 cm 8"
HI-04-590 23 cm 9"
 B/B Cvd.



$\frac{1}{4}$
HI-04-591 14.5 cm 5 3/4"
HI-04-592 18 cm 7"
HI-04-593 20 cm 8"
HI-04-594 23 cm 9"
 S/S Str.

Metzenbaum-Fino
Operating Scissors
T/C Insert



$\frac{1}{4}$
HI-04-595 14.5 cm 5 3/4"
HI-04-596 18 cm 7"
HI-04-597 20 cm 8"
HI-04-598 23 cm 9"
 S/S Cvd.



$\frac{1}{4}$
HI-04-599 14.5 cm 5 3/4"
HI-04-600 15.0 cm 6"
HI-04-601 17.0 cm 6 3/4"
HI-04-602 23.0 cm 9"
 B/B, Str.

Mayo Dissecting &
Gynecological Scissors
T/C Insert



$\frac{1}{4}$
HI-04-603 14.5 cm 5 3/4"
HI-04-604 15.0 cm 6"
HI-04-605 17.0 cm 6 3/4"
HI-04-606 23.0 cm 9"
 B/B, Cvd.



$\frac{1}{4}$
HI-04-607 15.0 cm 6"
HI-04-608 17.0 cm 6 3/4"
HI-04-609 19.5 cm 7 3/4"
HI-04-610 21.0 cm 8 3/4"
 B/B, Str.

Mayo-stille Dissecting &
Gynecological Scissors
T/C Insert



$\frac{1}{4}$
HI-04-611 15.0 cm 6"
HI-04-612 17.0 cm 6 3/4"
HI-04-613 19.5 cm 7 3/4"
HI-04-614 21.0 cm 8 3/4"
 B/B, Cvd.

OPERATING SCISSORS T.C.



Chirurgische scheren
Tijeras quirúrgicas

Ciseaux chirurgicaux
Forbici per chirurgia



ENT-DEAN SCISSORS
HI-04-615 17 cm, 6 3/4"
S/S Double Cvd.
Serrated



**KELLY OPERATING,
DISSECTING SCISSORS**

HI-04-616
16.0 cm, 6 1/4"
HI-04-617
18.0 cm, 7"
S/S Str.

HI-04-618
16.0 cm, 6 1/4"
HI-04-619
18.0 cm, 7"
S/S Cvd.



**Goldmann-Fox
Operating Scissors
T/C Insert**

1/4
HI-04-620
13 cm 5"
S/S Str.

1/4
HI-04-621
13 cm 5"
S/S Cvd.



WIRE CUTTING SCISSORS
TCHI-04-622 12 cm, 4 3/4"




**Lister Bandage Plaster
Scissors T/C Insert, one large ring**
HI-04-623 8.9cm, 3 1/2"
HI-04-624 11.43 cm, 4.5"
HI-04-625 14 cm, 5 1/2"
HI-04-626 16 cm, 6 1/4"
HI-04-627 18 cm, 7"
HI-04-628 20 cm, 8"



**Lister Bandage Plaster
Scissors, T/C Insert**
HI-04-629 8.9cm, 3 1/2"
HI-04-630 11.43 cm, 4.5"
HI-04-631 14 cm, 5 1/2"
HI-04-632 16 cm, 6 1/4"
HI-04-633 18 cm, 7"
HI-04-634 20 cm, 8"




HI-04-635
14 cm, 5 1/2"
B/B Str.



HI-04-636
14 cm, 5 1/2"
B/B Cvd.




HI-04-637
14 cm, 5 1/2"
S/B Str.



HI-04-638
14 cm, 5 1/2"
S/B Cvd.



HI-04-639
14 cm, 5 1/2"
S/S Str.



HI-04-640
14 cm, 5 1/2"
S/S Cvd.



HI-04-641
14 cm 5 1/4"
S/S Str.



HI-04-642
14 cm 5 1/4"
S/S Cvd.

**Castroviejo Scissors
T/C Insert**

**Joseph Operating
Scissors T/C Insert**



HI-04-643
16 cm 6 1/4"
HI-04-644
18 cm 7"
S/S Str.



HI-04-645
16 cm 6 1/4"
HI-04-646
18 cm 7"
S/S Cvd.



HI-04-648
10 cm 4"
HI-04-649
11 cm 4 1/4"
Str.




HI-04-650
10 cm 4"
HI-04-651
11 cm 4 1/4"
Cvd.



**Reynolds Operating
Scissors T/C Insert**

**Foman Operating
Scissors T/C Insert
HI-04-647
13 cm 5" Full Bent**

**Stevens Operating
Scissors T/C Insert**



HOMEDICA
INTERNATIONAL