


# Surgical Instruments Catalogue



HOMEDICA  
INTERNATIONAL



## **HOMEDICA INTERNATIONAL**

Homedica feels great honor introducing itself to be one of the leading manufacturing company of surgical & dental instruments. Professional always prefer Homedica instruments, which is our pride.

Our all instruments are manufactured by the high quality stainless steel made on our requirements, every instrument is hand crafted with the highest workmanship available in the industry. Through heat treatment, the instruments are hardened and tempered for proper functioning and long life. All instruments are polished electrolytically and manually for maximum resistance to corrosion and long life. Homedica pays special attention to the needs of doctors, with innovative designs to meet the needs of professionals. The wide assortment of standard patterns and exclusive Homedica designs gives you the most desirable range of surgical and dental instruments.

Over the several years, we have earned the reputation of for maintaining highest quality standards of quality and service. To manufacture instruments that meet with your changing demands is our achievement. Our highly skilled staff of professionals will continue to serve your needs and will try their outmost to uplift the standards which we have set in the past.

With your co-operation and guidance, Homedica will continue serving surgical and dental instruments in future.

**Management of Homedica international**

## BASIC METHODS FOR HANDLING AND CARE OF SURGICAL AND DENTAL INSTRUMENTS

Dear Customers:

We have received some questions concerning staining, pitting or rusting in stainless steel instruments. So, we took great honour and pleasure in guiding our respected customers to eliminate the problems about staining or rusting etc. The points are as follows:-

- 1) Rinse the instruments as soon as possible after use to remove blood and serum.
- 2) Dry the instruments after cleaning or before wrapping and storing, if not so, moisture left can cause rust and corrosion.
- 3) Do not leave instruments in cold sterilization solution for a long period, as the solution may lead the instrument to be corrosive.
- 4) Never use caustic cleaners, use those specified for cleaning surgical instruments.
- 5) Avoid cleaning instruments with steel wool; if want to clean with hand, use detergent soap and a soft brush and finally treat with distilled water.
- 6) You should check your autoclave; Iron, calcium, sodium or copper in water can cause rusting, staining or corrosion.
- 7) Leaving of instruments for a long time in the following solutions can also cause corrosion.
  - i) Phenol or Iodine.
  - ii) Aluminium, Barium or Calcium.
  - iii) Any Acid, Mercury or Potassium solution.
- 8) Try your utmost not to sterilize stainless steel instrument with chrome plated brass instruments, as any break in chrome surface may cause pitting or rusty surface in instruments.
- 9) Muslin wraps should be rinsed sufficiently by the laundry, so that all chlorine bleach may remove, which can cause corrosion or custing during autoclaving.

We hope that by acting on the above instructions, your problems relating to rusting, staining or pitting would abolish and the age of your instrument would increase.

**MANAGEMENT**

# HOLLOW WARES

MESH BASKETS

JARS

KIDNEY TRAYS

INSTRUMENTS BOXES

STERILIZING DRUMS





## STERILIZING DRUMS

	Ø	h
	mm	mm
HI-15-001	80	80
HI-15-002	120	120
HI-15-003	160	120
HI-15-004	190	140
HI-15-005	190	200
HI-15-006	240	120
HI-15-007	240	160
HI-15-008	240	240
HI-15-009	290	120
HI-15-010	290	160
HI-15-011	290	240
HI-15-012	340	160
HI-15-013	340	240
HI-15-014	390	145
HI-15-015	390	200
HI-15-016	390	290



## STERILIZING DRUM

	Ø	h
	mm	Size mm
HI-15-017	150	150
HI-15-018	175	175
HI-15-019	190	190
HI-15-020	200	140
HI-15-021	200	200
HI-15-022	240	120
HI-15-023	225	225
HI-15-024	240	240
HI-15-025	290	160
HI-15-026	280	240
HI-15-027	230	150
HI-15-028	300	250
HI-15-029	340	250
HI-15-030	350	250
HI-15-031	390	200



## STERILIZING DRUM

	Ø	h
	mm	Size mm
HI-15-032	150	150
HI-15-033	175	175
HI-15-034	190	190
HI-15-035	200	140
HI-15-036	200	200
HI-15-037	240	120
HI-15-038	225	225
HI-15-039	240	240
HI-15-040	290	160
HI-15-041	280	240
HI-15-042	230	150
HI-15-043	300	250
HI-15-044	340	250
HI-15-045	350	250
HI-15-046	390	200



## STERILIZING DRUM

	Ø	h
	mm	Size mm
HI-15-047	150	150
HI-15-048	175	175
HI-15-049	190	190
HI-15-050	200	140
HI-15-051	200	200
HI-15-052	240	120
HI-15-053	225	225
HI-15-054	240	240
HI-15-055	290	160
HI-15-056	280	240
HI-15-057	230	150
HI-15-058	300	250
HI-15-059	340	250
HI-15-060	350	250
HI-15-061	390	200



## STERILIZING DRUM

	Ø	h
	mm	Size mm
HI-15-062	150	150
HI-15-063	175	175
HI-15-064	190	190
HI-15-065	200	140
HI-15-066	200	200
HI-15-067	240	120
HI-15-068	225	225
HI-15-069	240	240
HI-15-070	290	160
HI-15-071	280	240
HI-15-072	230	150
HI-15-073	300	250
HI-15-074	340	250
HI-15-075	350	250
HI-15-076	390	200

# HOLLOWARE INSTRUMENTS



Hohlware Instrumenten | orfèverie de Instruments de  
Orfebrería Instrumentos | Holloware Gli strumenti



**G LOVES BOX**  
HI-15-077  
300 x 150 mm



**INSTRUMENT TRAY WITH COVER**  
HI-15-078



**INSTRUMENT TRAY**  
HI-15-079



**CATHETER TRAY**

	Size
HI-15-080	230 x 75 mm
HI-15-081	350 x 90 mm
HI-15-082	430 x 100 mm
HI-15-083	200 x 100 mm
HI-15-084	300 x 100 mm
HI-15-085	350 x 100 mm



**PERFORATED BOX**

	Size
HI-15-086	100 x 50 mm
HI-15-087	165 x 65 mm
HI-15-088	250 x 100 mm
HI-15-089	250 x 150 mm
HI-15-090	250 x 200 mm
HI-15-091	150 x 65 mm
HI-15-092	170 x 70 mm
HI-15-093	200 x 100 mm
HI-15-094	250 x 125 mm
HI-15-095	300 x 150 mm
HI-15-096	300 x 200 mm
HI-15-097	300 x 130 mm
HI-15-098	320 x 150 mm
HI-15-099	420 x 180 mm
HI-15-100	460 x 125 mm
HI-15-101	355 x 180 mm
HI-15-102	300 x 180 mm
HI-15-103	400 x 200 mm
HI-15-104	510 x 200 mm



**BOX WITH FILTER**

	Size
HI-15-105	100 x 50 mm
HI-15-106	165 x 65 mm
HI-15-107	250 x 100 mm
HI-15-108	250 x 150 mm
HI-15-109	250 x 200 mm
HI-15-110	150 x 65 mm
HI-15-111	170 x 70 mm
HI-15-112	200 x 100 mm
HI-15-113	250 x 125 mm
HI-15-114	300 x 150 mm
HI-15-115	300 x 200 mm
HI-15-116	300 x 130 mm
HI-15-117	320 x 150 mm
HI-15-118	420 x 180 mm
HI-15-119	460 x 125 mm
HI-15-120	355 x 180 mm
HI-15-121	300 x 180 mm
HI-15-122	400 x 200 mm
HI-15-123	510 x 200 mm




**INSTRUMENTS BOX**

Without Clip	Size
HI-15-124	100 x 50 mm
HI-15-125	165 x 65 mm
HI-15-126	250 x 100 mm
HI-15-127	250 x 150 mm
HI-15-128	250 x 200 mm
HI-15-129	150 x 65 mm
HI-15-130	170 x 70 mm
HI-15-131	200 x 100 mm
HI-15-132	250 x 125 mm
HI-15-133	300 x 150 mm
HI-15-134	300 x 200 mm
HI-15-135	300 x 130 mm
HI-15-136	320 x 150 mm
HI-15-137	420 x 180 mm
HI-15-138	460 x 125 mm
HI-15-139	355 x 180 mm
HI-15-140	300 x 180 mm
HI-15-141	400 x 200 mm
HI-15-142	510 x 200 mm



**INSTRUMENT BOX**

With Clip	Size
HI-15-143	100 x 50 mm
HI-15-144	165 x 65 mm
HI-15-145	250 x 100 mm
HI-15-146	250 x 150 mm
HI-15-147	250 x 200 mm
HI-15-148	150 x 65 mm
HI-15-149	170 x 70 mm
HI-15-150	200 x 100 mm
HI-15-151	250 x 125 mm
HI-15-152	300 x 150 mm
HI-15-153	300 x 200 mm
HI-15-154	300 x 130 mm
HI-15-155	320 x 150 mm
HI-15-156	420 x 180 mm
HI-15-157	460 x 125 mm
HI-15-158	355 x 180 mm
HI-15-159	300 x 180 mm
HI-15-160	400 x 200 mm
HI-15-161	510 x 200 mm



### INSTRUMENT BOILER

	Size
HI-15-162	180 x 80 x 40 mm
HI-15-163	200 x 100 x 50 mm
HI-15-164	250 x 125 x 40 mm
HI-15-165	300 x 150 x 60 mm
HI-15-166	350 x 175 x 75 mm
HI-15-167	400 x 200 x 100 mm



### PERFORATED TRAY

	Size
HI-15-168	180 x 80 x 40 mm
HI-15-169	200 x 100 x 50 mm
HI-15-170	250 x 125 x 40 mm
HI-15-171	300 x 150 x 60 mm
HI-15-172	350 x 175 x 75 mm
HI-15-173	400 x 200 x 100 mm



### INSTRUMENT BOILER WITH STAND

	Size
HI-15-174	180 x 80 x 40 mm
HI-15-175	200 x 100 x 50 mm
HI-15-176	250 x 125 x 40 mm
HI-15-177	300 x 150 x 60 mm
HI-15-178	350 x 175 x 75 mm
HI-15-179	400 x 200 x 100 mm



### INSTRUMENT BOILER ELECTRIC

	Size
HI-15-180	180 x 80 x 40 mm
HI-15-181	200 x 100 x 50 mm
HI-15-182	250 x 125 x 40 mm
HI-15-183	300 x 150 x 60 mm
HI-15-184	350 x 175 x 75 mm
HI-15-185	400 x 200 x 100 mm

# HOLLOWARE INSTRUMENTS




Hohlware Instrumenten  
Orfebreria Instrumentos

orfèverie de Instruments de  
Holloware Gli strumenti



**NEEDLE CASE**  
HI-15-186  
230 x 125 mm



**SCALER BOX**  
HI-15-187  
300 x 180 mm



**SCALER TRAY**  
HI-15-188  
100 x 200 mm



**INDO BOX WITH SILVER FRAME**  
HI-15-189



**Autoclaveable**  
HI-15-190



## INSTRUMENT TRAY WITHOUT COVER

	Size		Size
HI-15-191	240 x 240 mm	HI-15-202	230 x 125 mm
HI-15-192	355 x 400 mm	HI-15-203	200 x 150 mm
HI-15-193	300 x 150 mm	HI-15-204	250 x 200 mm
HI-15-194	300 x 180 mm	HI-15-205	300 x 250 mm
HI-15-195	300 x 240 mm	HI-15-206	350 x 300 mm
HI-15-196	340 x 320 mm	HI-15-207	400 x 350 mm
HI-15-197	230 x 180 mm	HI-15-208	450 x 300 mm
HI-15-198	250 x 130 mm	HI-15-209	300 x 200 mm
HI-15-199	200 x 130 mm	HI-15-210	350 x 250 mm
HI-15-200	180 x 130 mm	HI-15-211	250 x 125 mm
HI-15-201	150 x 100 mm	HI-15-212	200 x 100 mm
		HI-15-213	400 x 200 mm

## INSTRUMENT TRAY WITHOUT COVER

	Size		Size
HI-15-214	240 x 240 mm	HI-15-225	230 x 125 mm
HI-15-215	355 x 400 mm	HI-15-226	200 x 150 mm
HI-15-216	300 x 150 mm	HI-15-227	250 x 200 mm
HI-15-217	300 x 180 mm	HI-15-228	300 x 250 mm
HI-15-218	300 x 240 mm	HI-15-229	350 x 300 mm
HI-15-219	340 x 320 mm	HI-15-230	400 x 350 mm
HI-15-220	230 x 180 mm	HI-15-231	450 x 300 mm
HI-15-221	250 x 130 mm	HI-15-232	300 x 200 mm
HI-15-222	200 x 130 mm	HI-15-233	350 x 250 mm
HI-15-223	180 x 130 mm	HI-15-234	250 x 125 mm
HI-15-224	150 x 100 mm	HI-15-235	200 x 100 mm
		HI-15-236	400 x 200 mm





**KIDNEY TRAY WITH COVER**  
(English Pattern)

	Size
HI-15-237	150 mm
HI-15-238	200 mm
HI-15-239	240 mm
HI-15-240	300 mm



**KIDNEY TRAY (German Pattern)**  
(without cover)

	Size
HI-15-241	150 mm
HI-15-242	200 mm
HI-15-243	240 mm
HI-15-244	300 mm



**SPLIT LAMP**  
HI-15-245  
100 ml.



**BOWL**

	Size
HI-15-246	80 mm
HI-15-247	100 mm
HI-15-248	125 mm
HI-15-249	150 mm
HI-15-250	200 mm
HI-15-251	250 mm
HI-15-252	300 mm
HI-15-253	350 mm
HI-15-254	400 mm



**WASH BASIN**

	Size
HI-15-255	80 mm
HI-15-256	100 mm
HI-15-257	125 mm
HI-15-258	150 mm
HI-15-259	200 mm
HI-15-260	250 mm
HI-15-261	300 mm
HI-15-262	350 mm
HI-15-263	400 mm



**BOWL**

	Size
HI-15-264	80 mm
HI-15-265	100 mm
HI-15-266	125 mm
HI-15-267	150 mm
HI-15-268	200 mm
HI-15-269	250 mm
HI-15-270	300 mm
HI-15-271	350 mm
HI-15-272	400 mm



**DRESSING JAR WITH KNOB**

	Size
HI-15-273	80 x 80 mm
HI-15-274	100 x 100 mm
HI-15-275	125 x 125 mm
HI-15-276	150 x 150 mm
HI-15-277	200 x 200 mm



**DRESSING JAR**

	Size
HI-15-278	80 x 80 mm
HI-15-279	100 x 100 mm
HI-15-280	125 x 125 mm
HI-15-281	150 x 150 mm
HI-15-282	200 x 200 mm



**DRESSING JAR**

	Size
HI-15-283	80 x 80 mm
HI-15-284	100 x 100 mm
HI-15-285	125 x 125 mm
HI-15-286	150 x 150 mm
HI-15-287	200 x 200 mm



**DRESSING JAR**

	Size
HI-15-288	80 x 80 mm
HI-15-289	100 x 100 mm
HI-15-290	125 x 125 mm
HI-15-291	150 x 150 mm
HI-15-292	200 x 200 mm

# HOLLOWARE INSTRUMENTS



Hohlware Instrumenten  
Orfebreria Instrumentos

orfèverie de Instruments de  
Holloware Gli strumenti



**FORCEPS  
JAR**

	Size
HI-15-293	55 x 100 mm
HI-15-294	55 x 140 mm
HI-15-295	55 x 170 mm
HI-15-296	55 x 180 mm
HI-15-297	80 x 100 mm



**FORCEPS  
JAR**

	Size
HI-15-298	55 x 100 mm
HI-15-299	55 x 140 mm
HI-15-300	55 x 170 mm
HI-15-301	55 x 180 mm
HI-15-302	80 x 100 mm



**FORCEPS  
JAR**

	Size
HI-15-303	55 x 100 mm
HI-15-304	55 x 140 mm
HI-15-305	55 x 170 mm
HI-15-306	55 x 180 mm
HI-15-307	80 x 100 mm



**FORCEPS  
JAR**

	Size
HI-15-308	55 x 100 mm
HI-15-309	55 x 140 mm
HI-15-310	55 x 170 mm
HI-15-311	55 x 180 mm
HI-15-312	80 x 100 mm



**FORCEPS  
JAR**

	Size
HI-15-313	55 x 100 mm
HI-15-314	55 x 140 mm
HI-15-315	55 x 170 mm
HI-15-316	55 x 180 mm
HI-15-317	80 x 100 mm



**COTTON HOLDER**  
HI-15-318



**GALIPOTS**

	Capacity
HI-15-319	2 Ounces
HI-15-320	3 Ounces
HI-15-321	4 Ounces
HI-15-322	5 Ounces
HI-15-323	6 Ounces
HI-15-324	7 Ounces
HI-15-325	8 Ounces



**GALIPOTS**

	Capacity
HI-15-326	2 Ounces
HI-15-327	3 Ounces
HI-15-328	4 Ounces
HI-15-329	5 Ounces
HI-15-330	6 Ounces
HI-15-331	7 Ounces
HI-15-332	8 Ounces



**SPRIT LAMP**  
HI-15-333  
100 ml.



**INDO BOX ROUND**  
HI-15-334



**NEEDLE CASE**  
HI-15-335



**MEASURING JUG**  
Capacity  
HI-15-336 1.00 Ltr  
HI-15-337 1.50 Ltr  
HI-15-338 2.00 Ltr



**MEDICINE CUP GRADUATED**  
Capacity  
HI-15-339 50 ml  
HI-15-340 60 ml  
HI-15-341 100 ml



**BUCKET WITH LID**  
Capacity  
HI-15-342 10 Ltr  
HI-15-343 12 Ltr  
HI-15-344 15 Ltr



**SPUTUM MUG**  
**HI-15-345**  
 100 x 50 mm



**ENEMA CANE**  
 Capacity  
**HI-15-346** 0.5 Ltr  
**HI-15-347** 1.5 Ltr  
**HI-15-348** 2.0 Ltr



**DUST BIN**  
**HI-15-349**  
 20 Lit



**ROUND BEDPAN**  
**HI-15-350**  
 300 x 85 mm



**BEDPAN**  
**HI-15-351**  
 Size: Adult, Child



**BEDPAN**  
**HI-15-352**  
 Size: Adult, Child



**FUNNAL**  
**HI-15-353**  
 100 mm



**URINAL - FEMALE**  
**HI-15-354**




**URINAL - MALE**  
**HI-15-355**



**SHAVING MUG**  
**HI-15-356**  
 100 mm


## PERFORATED BASKETS

perforiert Körbe | paniers perforés  
Cestas perforadas | ceste forate



### MESH BASKET

HI-15-357 240 x 258 x 50 mm  
HI-15-358 250 x 240 x 100 mm  
HI-15-359 240 x 250 x 60 mm



### MESH BASKET


HI-15-360 480 x 258 x 50 mm  
HI-15-361 480 x 258 x 100 mm  
HI-15-362 520 x 258 x 60 mm  
HI-15-363 520 x 258 x 70 mm  
HI-15-364 450 x 340 x 70 mm  
HI-15-365 320 x 230 x 70 mm



### MESH BASKET


With Perforated Sheet Side

HI-15-366 230 x 130 x 50 mm  
HI-15-367 230 x 130 x 60 mm  
HI-15-368 240 x 130 x 50 mm  
HI-15-369 240 x 130 x 60 mm  
HI-15-370 240 x 250 x 50 mm  
HI-15-371 240 x 250 x 60 mm  
HI-15-372 480 x 250 x 50 mm  
HI-15-373 480 x 250 x 60 mm  
HI-15-374 480 x 250 x 70 mm  
HI-15-375 300 x 300 x 50 mm  
HI-15-376 350 x 350 x 60 mm  
HI-15-377 470 x 400 x 60 mm  
HI-15-378 480 x 480 x 50 mm  
HI-15-379 480 x 480 x 60 mm  
HI-15-380 500 x 300 x 50 mm  
HI-15-381 520 x 300 x 50 mm  
HI-15-382 520 x 300 x 60 mm  
HI-15-383 550 x 430 x 50 mm  
HI-15-384 600 x 300 x 60 mm



### ARTHROSCOPE BASKET

HI-15-385 290 x 80 x 52 mm



### ENDOSCOPE BASKET

HI-15-386 460 x 80 x 52 mm  
HI-15-387 660 x 80 x 52 mm  
HI-15-388 670 x 80 x 52 mm


## STERILIZING FORCEPS

Sterilisieren Zange  
Pinzas de esterilización

Pinces stérilisation  
Pinze di sterilizzazione



**DAVIS**  
**HI-15-389** 16 cm, 6 ¼"  
**HI-15-390** 20 cm, 8"  
**HI-15-391** 24 cm, 9 ½"



**UTILITY**  
**HI-15-392** 20 cm, 8"  
**HI-15-393** 28 cm, 11"



**HI-15-394** 20 cm, 8"  
**HI-15-395** 25 cm, 10"  
**HI-15-396** 30 cm, 11 ¾"




**ROGGE**  
**HI-15-397**  
 27 cm, 11"



**CHEATEL**  
 Box Joint  
**HI-15-398**  
 27 cm, 11"



**CHEATEL**  
 Screw Joint  
**HI-15-399**  
 27 cm, 11"



**HOMEDICA**  
INTERNATIONAL