


Surgical Instruments Catalogue



HOMEDICA
INTERNATIONAL



HOMEDICA INTERNATIONAL

Homedica feels great honor introducing itself to be one of the leading manufacturing company of surgical & dental instruments. Professional always prefer Homedica instruments, which is our pride.

Our all instruments are manufactured by the high quality stainless steel made on our requirements, every instrument is hand crafted with the highest workmanship available in the industry. Through heat treatment, the instruments are hardened and tempered for proper functioning and long life. All instruments are polished electrolytically and manually for maximum resistance to corrosion and long life. Homedica pays special attention to the needs of doctors, with innovative designs to meet the needs of professionals. The wide assortment of standard patterns and exclusive Homedica designs gives you the most desirable range of surgical and dental instruments.

Over the several years, we have earned the reputation of for maintaining highest quality standards of quality and service. To manufacture instruments that meet with your changing demands is our achievement. Our highly skilled staff of professionals will continue to serve your needs and will try their outmost to uplift the standards which we have set in the past.

With your co-operation and guidance, Homedica will continue serving surgical and dental instruments in future.

Management of Homedica international

BASIC METHODS FOR HANDLING AND CARE OF SURGICAL AND DENTAL INSTRUMENTS

Dear Customers:

We have received some questions concerning staining, pitting or rusting in stainless steel instruments. So, we took great honour and pleasure in guiding our respected customers to eliminate the problems about staining or rusting etc. The points are as follows:-

- 1) Rinse the instruments as soon as possible after use to remove blood and serum.
- 2) Dry the instruments after cleaning or before wrapping and storing, if not so, moisture left can cause rust and corrosion.
- 3) Do not leave instruments in cold sterilization solution for a long period, as the solution may lead the instrument to be corrosive.
- 4) Never use caustic cleaners, use those specified for cleaning surgical instruments.
- 5) Avoid cleaning instruments with steel wool; if want to clean with hand, use detergent soap and a soft brush and finally treat with distilled water.
- 6) You should check your autoclave; Iron, calcium, sodium or copper in water can cause rusting, staining or corrosion.
- 7) Leaving of instruments for a long time in the following solutions can also cause corrosion.
 - i) Phenol or Iodine.
 - ii) Aluminium, Barium or Calcium.
 - iii) Any Acid, Mercury or Potassium solution.
- 8) Try your utmost not to sterilize stainless steel instrument with chrome plated brass instruments, as any break in chrome surface may cause pitting or rusty surface in instruments.
- 9) Muslin wraps should be rinsed sufficiently by the laundry, so that all chlorine bleach may remove, which can cause corrosion or custing during autoclaving.

We hope that by acting on the above instructions, your problems relating to rusting, staining or pitting would abolish and the age of your instrument would increase.

MANAGEMENT

SCALPEL HANDLE, SCALPEL BLADES, DISSECTING KNIVES & SAWS



SCALPEL HANDLES



Skalpеллgriffe
Mangos de bisturias

Manches pour bistouris
Manici per bisturi



Locks At 90°, 60°, 30°,
Straight, 120°, 150°, 180°,
Angles



SCALPEL HANDLE #3
HI-03-001
12.75 cm, 5"
Blade Size 10, 11,
12, 13, 15, 17



**#3 SCALPEL HANDLE/
GRADUATED**
HI-03-002
12.75 cm, 5"
Blade Size 10, 11,
12, 13, 15, 17



**SCALPEL HANDLE
#4 GRADUATED**
HI-03-003
with scale



SCALPEL HANDLE #4
HI-03-004
12.75 cm, 5" Blade
Size 18, 20, 21, 22,
23, 24, 25, 60,70



SCALPEL HANDLE 4LA
HI-03-005
21 cm, 8 1/4" Blade
Size 18, 20, 21, 22,
23, 24, 25, 60,70



**SIEGEL ROUND
SCALPEL HANDLE
#3**
HI-03-006
With Adjustable
Tip 16 cm, 6 1/4"
Blade Size 10,
11, 12, 13, 15,
17

SCALPEL HANDLES

Skalpgriffe
Mangos de bisturíes

Manches pour bistouris
Manici per bisturi



SCALPEL HANDLE #3, #4
HI-03-007
Double Ended
16 cm, 6"

SCALPEL HANDLE 3L LONG
HI-03-008
21 cm 8 1/4"
Blade Size
10, 11, 12, 15, 17

SCALPEL HANDLE #7
HI-03-009
16.5 cm, 6 1/4"
Blade 10, 11,
12, 13, 15, 17

SCALPEL HANDLE 4L LONG
HI-03-010
21 cm, 8 1/4"
Blade
Size 18, 20, 21,
22, 23, 24, 25,
60,70

SCALPEL HANDLE / LONG ANGLED #3LA
HI-03-011
20.25 cm, 8 1/4"
Blade Size
10, 11, 12, 15, 17

KNIFE HANDLE
HI-03-012
(Smooth) for
Blades Nos.
60 and 70

ROUND KNURLED SCALPEL HANDLE
HI-03-013
12.7 cm, 5"
Blade Size
50, 60 & 70

SCALPEL HANDLES

Skalpelligriffe
Mangos de bisturios

Manches pour bistouris
Manici per bisturi



**ROUND KNURLED
SCALPEL HANDLE**
HI-03-014
12 cm, 4 3/4"



**KAYE ROUND
SCALPEL HANDLE**
#3RK
HI-03-015
15 cm, 5 3/4", Blade
Size 10, 11,
12, 13, 15, 17



**CORWIN ROUND
SCALPEL HANDLE**
#3R
HI-03-016
14.5 cm, 5 1/2"
Blade Size 10,
11, 12, 15, 17



SCALPEL HANDLE
#3RC
HI-03-017
14 cm, 5 1/2"
Blade Size 10,
11, 12, 13, 15,
17



**CORWIN ROUND
SCALPEL HANDLE,**
#4R
HI-03-018
14.5 cm, 5 1/2",
Blade Size 18, 20,
21, 22, 23, 24,
25, 60, 70



SCALPEL HANDLE
#3P
HI-03-019
15.5 cm,
6", Blade Size
10, 11, 12, 13,
15, 17

SCALPEL BLADES

Skalpellskliven, | Hojas de bisturi, Manches
Lames de bistouris, Mangos | Lame per bisturi, Manici



HI-03-020
Sterile Scalpel
Blade Size # 6
Stainless Steel



HI-03-021
Sterile Scalpel
Blade Size # 10
Stainless Steel



HI-03-022
Sterile Scalpel
Blade Size # 10A
Stainless Steel



HI-03-023
Sterile Scalpel
Blade Size # 11
Stainless Steel



HI-03-024
Sterile Scalpel
Blade Size # 11P
Stainless Steel



HI-03-025
Sterile Scalpel
Blade Size # E11
Stainless Steel



HI-03-026
Sterile Scalpel
Blade Size # 12
Stainless Steel



HI-03-027
Sterile Scalpel
Blade Size # 13
Stainless Steel



HI-03-028
Sterile Scalpel
Blade Size # 15
Stainless Steel



HI-03-029
Sterile Scalpel
Blade Size # 15A
Stainless Steel



HI-03-030
Sterile Scalpel
Blade Size # 17
Stainless Steel



HI-03-031
Sterile Scalpel
Blade Size # 18
Stainless Steel



HI-03-032
Sterile Scalpel
Blade Size # 20
Stainless Steel



HI-03-033
Sterile Scalpel
Blade Size # 21
Stainless Steel



HI-03-034
Sterile Scalpel
Blade Size # 22
Stainless Steel



HI-03-035
Sterile Scalpel
Blade Size # 23
Stainless Steel



HI-03-036
Sterile Scalpel
Blade Size # 24
Stainless Steel



HI-03-037
Sterile Scalpel
Blade Size # 25
Stainless Steel



HI-03-038
Sterile Scalpel
Blade Size # 36
Stainless Steel



HI-03-039
Sterile Scalpel
Blade Size # 40
Stainless Steel



HI-03-040
Sterile Scalpel
Blade Size # 60
Stainless Steel



HI-03-041
Sterile Scalpel
Blade Size # 70
Stainless Steel

SCALPEL HANDLES



Skalpelligriffe
Mangos de bisturios

Manches pour bistouris
Manici per bisturi



DIEFFENBACH
16 cm, 6 1/4"



HI-03-042
Fig. 1
16 cm, 6 1/4"



HI-03-043
Fig. 2
16 cm, 6 1/4"



HI-03-044
Fig. 3
16 cm, 6 1/4"



HI-03-045
Fig. 4
16 cm, 6 1/4"



HI-03-046
Fig. 5
16 cm, 6 1/4"



HI-03-047
Fig. 6
16 cm, 6 1/4"



HI-03-048
Fig. 7
16 cm, 6 1/4"



HI-03-049
Fig. 8
16 cm, 6 1/4"



HI-03-050
Fig. 9
16 cm, 6 1/4"



HI-03-051
Fig. 10
16 cm, 6 1/4"



HI-03-052
Fig. 15
16 cm, 6 1/4"




HI-03-053
Fig. 16
16 cm, 6 1/4"




COLLIN DISSECTING KNIVES Blade Size
 HI-03-054 13 cm, 5"
 HI-03-055 16 cm, 6 ¼"
 HI-03-056 19 cm, 7 ½"
 HI-03-057 22 cm, 8 ½"
 HI-03-058 28 cm, 11 ¼"



CATLIN DISSECTING KNIVES Blade Size
 HI-03-059 16 cm, 6 ¼"
 HI-03-060 19 cm, 7 ½"
 HI-03-061 22 cm, 8 ½"



LISTON DISSECTING KNIVES Blade Size
 HI-03-062 16 cm, 6 ¼"
 HI-03-063 19 cm, 7 ½"
 HI-03-064 22 cm, 8 ½"



LANGENBECK DISSECTING KNIVES BLADE
 HI-03-065
 12 cm, 4 ¾"



DIEFFENBACH DISSECTING KNIVES
 HI-03-066
 13 cm, 5"



DIEFFENBACH DISSECTING KNIVES
 HI-03-067
 13 cm, 5"



SIMON DISSECTING KNIVES
 23.5 cm, 9 ½"



HI-03-068
Fig. 1



HI-03-069
Fig. 2



HI-03-070
Fig. 3



HI-03-071
Fig. 4

AMPUTATION KNIVES & RESECTION SAWS



Amputationsmesser, Cuchillos para amputación
Couteaux à amputation, coltelli per amputazione

Resektionsmesser, Cuchillos para resección
Couteaux à résection, Coltelli per resezione



AMPUTATION KNIFE
HI-03-072
29 cm, 11 1/2"



AMPUTATION KNIFE
HI-03-073
32.5 cm, 13"



CATLIN AMPUTATION
HI-03-074
29 cm, 11 1/2"



RESECTION KNIFE
HI-03-075
18 cm, 7"




LANGENBECK AMPUTATION
HI-03-076
18 cm, 7"



LANGENBECK AMPUTATION
HI-03-077
18 cm, 7"



RAZOR HOLLOW
Blades On
Both Sides
Blade Flat / Hollow
HI-03-078



CHARRIERE AMPUTATION & RESECTION
HI-03-080 27 cm, 11"
HI-03-081 30 cm, 12"



CHARRIERE AMPUTATION & RESECTION
HI-03-082
35 cm, 13 3/4"



SATTERLEE AMPUTATION & RESECTION
HI-03-083
31 cm, 12 1/4"



HI-03-085
spare blade


HUMBY DERMATOMES
HI-03-084 31cm 12 1/4"
Complete with 1 Blade



WATSON DERMATOMES
HI-03-086 30cm 12"
Complete with 1 Blade,
fixed bar



SILVER DERMATOMES
HI-03-087 19cm 7 1/2"



HOMEDICA
INTERNATIONAL